


# A Robust Toolkit for Functional Profiling of the Yeast Genome



Xuwen Pan, Daniel S. Yuan, Dong Xiang, Xiaoling Wang, Sharon Sookhai-Mahadeo, Joel S. Bader, Philip Hieter, Forrest Spencer, and Jef D. Boeke



**Figure S1. The SGA Reporter Can Be Exploited to Convert Diploid into Haploid Yeast**

(A) An experimental design to test the efficiency of haploid selection with the SGA reporter. A heterozygous diploid YKO mutant of an essential gene (*ESSX/essxΔ::kanMX*) containing the SGA reporter is sporulated. Haploid progenies are selected on the magic medium with (MM) or without (MM lacking G418) G418. Cells grown on the MM medium should contain both the wild-type and the *kanMX* alleles of the essential gene. These *kanMX* allele-containing cells likely result from meiotic nondisjunction or more rarely from mitotic homozygosity of both the mating-type and *CAN1* loci. These are not true haploids and we call them pseudo-haploid cells. By comparing the number of wild-type haploid cells and the pseudo-haploid cells generated from the same sporulation one can estimate the power of haploid selection with the SGA reporter.


(B) The haploid-selection power of the SGA reporter is around  $10^3$ . Three "SGA reporter"-containing heterozygous diploid YKO mutants of essential genes (*CDC102/cdc102Δ::kanMX*, *XPY452a/α*; *RPF2/rpf2Δ::kanMX*, *XPY454a/α*; *YKR079C/ykr079cΔ::kanMX*, *XPY460a/α*) were sporulated. Meiotic progenies were selected on the magic medium with or without G418. The number of colonies from the same amount of sporulation culture generated on both media was compared. A ratio of true haploid (without G418) over pseudo-haploid (with G418) was plotted.



**Figure S2.** High Accuracy of Integrative Transformation in the Heterozygous Diploid YKOs

(A) A diagram of the test of transformation accuracy in a heterozygous diploid YKO mutant. A gene disruption cassette of an essential locus (*essxΔ::natMX* with ~1.5kb flanking sequence) was transformed into a preexisting heterozygous YKO diploid mutant of the same essential gene (*ESSX/essxΔ::kanMX*). When the *essxΔ::natMX* cassette integrates into the right locus, it either disrupt the wild-type *ESSX* gene and cause lethality or replace the *essxΔ::kanMX* copy and the resultant transformants become resistant to CloNat but sensitive to G418 (*Nat<sup>R</sup> G418<sup>s</sup>*). When ectopic integration happens, the resultant transformants are resistant to both drugs (*Nat<sup>R</sup> G418<sup>R</sup>*). By comparing the number of *Nat<sup>R</sup> G418<sup>s</sup>* and *Nat<sup>R</sup> G418<sup>R</sup>* colonies from the same transformation reaction, one can calculate what portion of the *Nat<sup>R</sup>* transformants are results of accurate integration of the knockout cassette into that essential locus in the heterozygous diploid strain.

(B) Integrative transformation accuracy in heterozygous diploid YKOs. Four different gene-disruption cassettes depicted were individually transformed into the preexisting heterozygous diploid YKO mutants of the corresponding essential genes. Transformation accuracy (in the form of percentage) for each reaction was plotted.



**Figure S3. The Heterozygous Diploid Pool Is Easy to Transform**

(A) A diagram of making a query construct for the integrative transformation. A *URA3* marker was used to replace the pre-existing *xxxΔ::kanMX* module in each YKO mutant (haploid or diploid) via homologous recombination (*xxx* represents any YKO mutation). The resultant *Ura*<sup>+</sup> G418<sup>r</sup> mutant is used as the template in a long-range PCR reaction to amplify the query construct (*xxxΔ::URA3*) for integrative transformation. An *xxxΔ::natMX* cassette can be constructed in a similar manner.

(B) High integrative transformation yields in a heterozygous diploid pool. Twelve randomly chosen gene-disruption cassettes depicted were individually en masse transformed into a pool of 5,896 heterozygous diploid YKO mutants containing the SGA reporter as described in Materials and Methods. The number of *Ura*<sup>+</sup> transformants obtained from each reaction was plotted.

**Table S1. Behaviours of DNA damage response YKO mutants in high throughput screens.**

| Strain<br>genotype | Type of YKOs used                                     |          |          |                                      |          |    |  |    |    |
|--------------------|---|----------|----------|--------------------------------------|----------|----|--|----|----|
|                    | <i>MATa</i> /□ homozygous<br>diploid YKO <sup>a</sup> |          |          | <i>MATa</i> haploid YKO <sup>b</sup> |          |    | heterozygous diploid<br>YKO derivatives <sup>c</sup> |    |    |
|                    | MMS   | HU       | UV       | MMS                                  | HU       | UV | MMS  | HU | UV |
| <i>ddc1</i> □      | NA  | NA       | NA       | S                                    | S        | S  | S  | S  | S  |
| <i>mec3</i> □      | <b>R</b>  | <b>R</b> | S        | S                                    | S        | S  | S  | S  | S  |
| <i>rad17</i> □     | S   | <b>R</b> | S        | S                                    | S        | S  | S  | S  | S  |
| <i>rad24</i> □     | NA  | NA       | NA       | S                                    | S        | S  | S  | S  | S  |
| <i>rad9</i> □      | <b>R</b>  | <b>R</b> | <b>R</b> | S                                    | <b>R</b> | S  | S  | S  | S  |

Notes: **R**- resistant; S-sensitive; NA-no data is available. Data in this table was adopted from: <sup>a</sup> Bennett et al. (2001); <sup>b</sup> Chang et al., (2002); <sup>c</sup> Pan et al., in preparation.

**Table S2. Benomyl-sensitivity of YKO mutants studied by microarray analysis**

| Gene Name | Benomyl concentration ( $\mu$ g/ml) |       |        |        |        |        |        |        |        |
|-----------|-------------------------------------|-------|--------|--------|--------|--------|--------|--------|--------|
|           | 1                                   | 5     | 10     | 15     | 20     | 25     | 30     | 35     | 40     |
| CIN1      | 4.28                                | 83.68 | 100.61 | 148.13 | 98.20  | 371.48 | 176.94 | 212.63 | 169.98 |
| YML094C-A | 2.66                                | 18.12 | 47.09  | 30.68  | 26.18  | 77.70  | 35.66  | 28.90  | 30.00  |
| PAC10     | 2.63                                | 9.28  | 51.57  | 75.01  | 19.63  | 75.57  | 19.73  | 7.28   | 32.89  |
| PFD1      | 2.59                                | 14.63 | 60.19  | 71.96  | 33.03  | 101.15 | 110.11 | 57.04  | 30.98  |
| GIM3      | 2.37                                | 7.86  | 47.91  | 84.08  | 26.35  | 94.88  | 33.27  | 39.20  | 29.48  |
| TUB3      | 2.28                                | 23.01 | 80.12  | 79.25  | 93.04  | 110.69 | 87.58  | 77.60  | 107.84 |
| GIM5      | 2.14                                | 17.88 | 44.67  | 28.74  | 17.25  | 118.62 | 21.87  | 50.97  | 29.16  |
| YKE2      | 1.74                                | 8.36  | 28.07  | 29.29  | 7.56   | 29.18  | 26.05  | 25.05  | 25.05  |
| GIM4      | 1.58                                | 5.55  | 25.76  | 85.44  | 319.30 | 407.87 | 350.77 | 278.47 | 421.53 |
| CIN2      | 1.49                                | 74.87 | 38.18  | 56.57  | 72.59  | 147.35 | 121.31 | 77.04  | 54.52  |
| NFU1      | 1.15                                | 5.37  | 5.09   | 8.72   | 17.77  | 41.56  | 53.82  | 31.13  | 85.97  |
| WHI2      | 1.10                                | 6.77  | 10.23  | 2.66   | 4.32   | 67.88  | 41.75  | 183.34 | 60.69  |
| SOD2      | 1.00                                | 5.84  | 32.29  | 19.03  | 28.67  | 100.02 | 47.70  | 23.45  | 264.13 |
| ISU1      | 0.97                                | 8.84  | 9.65   | 11.55  | 13.95  | 102.81 | 57.09  | 43.87  | 10.39  |
| CIN4      | 0.86                                | 0.82  | 5.27   | 21.26  | 73.96  | 83.90  | 47.97  | 56.54  | 129.72 |
| CCR4      | 0.84                                | 0.90  | 1.85   | 1.14   | 2.47   | 14.96  | 18.94  | 138.82 | 116.63 |
| PYC1      | 0.89                                | 6.00  | 9.64   | 16.73  | 15.94  | 13.29  | 39.59  | 31.28  | 36.72  |
| UBA3      | 1.09                                | 4.21  | 7.89   | 10.10  | 21.18  | 31.58  | 27.70  | 13.66  | 14.86  |
| POP2      | 0.83                                | 0.97  | 1.26   | 1.97   | 6.38   | 35.36  | 36.64  | 31.98  | 95.14  |
| MIR1      | 1.15                                | 11.36 | 47.08  | 9.26   | 6.54   | 9.33   | 9.70   | 13.93  | 12.96  |
| UBI4      | 1.36                                | 0.28  | 0.91   | 7.86   | 4.65   | 32.59  | 16.18  | 25.94  | 24.87  |
| DHH1      | 1.13                                | 0.53  | 1.10   | 1.22   | 3.59   | 12.79  | 34.05  | 30.00  | 19.85  |
| YBR030W   | 0.75                                | 1.40  | 1.14   | 1.02   | 1.09   | 3.02   | 4.70   | 67.94  | 232.45 |
| THR1      | 0.81                                | 1.95  | 13.33  | 17.75  | 15.90  | 6.08   | 10.17  | 15.00  | 20.31  |
| RTG1      | 0.72                                | 2.80  | 7.55   | 7.10   | 13.46  | 9.03   | 10.64  | 28.27  | 21.12  |
| BRO1      | 1.47                                | 1.37  | 1.65   | 5.59   | 14.19  | 34.46  | 10.14  | 8.47   | 25.80  |
| RTG2      | 1.02                                | 3.15  | 9.62   | 9.22   | 14.19  | 9.21   | 12.14  | 14.16  | 22.75  |
| MAD2      | 1.20                                | 1.30  | 1.26   | 1.04   | 2.36   | 6.86   | 21.20  | 36.79  | 227.73 |
| YJR018W   | 1.05                                | 2.80  | 2.58   | 1.71   | 3.98   | 2.96   | 6.50   | 49.56  | 35.45  |
| YGR046W   | 0.86                                | 6.37  | 3.04   | 3.28   | 9.30   | 19.82  | 8.09   | 19.72  | 42.06  |
| FYV6      | 1.43                                | 1.97  | 7.85   | 11.27  | 7.66   | 16.85  | 9.87   | 13.56  | 10.19  |
| PPZ1      | 1.03                                | 4.32  | 4.14   | 2.77   | 5.53   | 12.49  | 13.07  | 26.33  | 28.36  |
| YDR334W   | 1.19                                | 1.24  | 1.11   | 2.24   | 4.46   | 8.57   | 12.84  | 34.67  | 24.13  |
| SNT309    | 0.89                                | 1.92  | 0.86   | 1.32   | 2.41   | 14.49  | 18.48  | 25.79  | 11.82  |
| YOR073W   | 1.22                                | 2.64  | 3.80   | 5.19   | 8.17   | 13.32  | 6.31   | 25.39  | 11.93  |
| PSO2      | 0.83                                | 0.92  | 1.34   | 1.15   | 1.43   | 9.09   | 24.09  | 27.05  | 76.23  |
| RMD11     | 0.95                                | 1.97  | 2.46   | 1.59   | 2.39   | 7.79   | 16.67  | 28.17  | 60.22  |
| AFG3      | 1.09                                | 3.42  | 3.03   | 1.84   | 3.13   | 6.05   | 8.76   | 33.99  | 15.89  |
| GPB2      | 1.18                                | 1.40  | 1.26   | 1.32   | 1.14   | 23.95  | 15.15  | 14.61  | 28.27  |
| GLR1      | 1.05                                | 3.80  | 3.45   | 3.62   | 5.15   | 20.70  | 8.31   | 12.46  | 20.80  |

|         |      |       |       |      |       |       |       |       |       |
|---------|------|-------|-------|------|-------|-------|-------|-------|-------|
| YAL046C | 1.19 | 2.46  | 4.83  | 6.24 | 5.78  | 13.65 | 10.19 | 12.20 | 20.57 |
| CTK1    | 1.38 | 0.67  | 1.15  | 1.86 | 5.85  | 16.46 | 12.81 | 16.12 | 22.58 |
| YKL134C | 0.87 | 5.25  | 3.27  | 2.27 | 3.61  | 6.20  | 16.49 | 18.30 | 29.74 |
| GPB1    | 1.05 | 1.83  | 1.25  | 0.84 | 1.13  | 24.01 | 8.99  | 16.69 | 17.32 |
| SKN7    | 0.84 | 1.08  | 2.12  | 6.86 | 17.41 | 15.85 | 6.38  | 5.21  | 3.99  |
| RTG3    | 1.09 | 2.32  | 12.09 | 9.17 | 7.83  | 6.96  | 8.34  | 6.99  | 21.53 |
| YLR294C | 1.09 | 4.87  | 4.16  | 1.85 | 3.16  | 4.55  | 7.70  | 26.79 | 14.16 |
| EAP1    | 1.21 | 1.48  | 1.21  | 1.16 | 2.86  | 7.95  | 10.25 | 27.80 | 26.53 |
| YME1    | 1.07 | 4.44  | 1.31  | 1.17 | 2.82  | 9.80  | 11.32 | 19.83 | 16.25 |
| RIM15   | 0.92 | 0.50  | 0.73  | 0.75 | 0.75  | 14.62 | 18.64 | 13.82 | 10.16 |
| POR1    | 1.01 | 3.13  | 1.99  | 1.45 | 2.41  | 4.69  | 7.81  | 26.57 | 29.97 |
| DBF2    | 1.13 | 11.66 | 2.90  | 0.87 | 3.53  | 4.85  | 6.89  | 16.82 | 23.98 |
| PAN2    | 1.08 | 1.16  | 1.73  | 2.18 | 3.49  | 4.03  | 11.10 | 23.20 | 71.31 |
| MCT1    | 0.65 | 3.70  | 2.51  | 1.20 | 2.14  | 6.10  | 15.72 | 15.10 | 12.46 |
| CTK3    | 1.67 | 4.12  | 2.38  | 2.78 | 7.32  | 11.42 | 7.18  | 9.40  | 21.29 |
| RPN4    | 0.81 | 1.69  | 1.63  | 4.28 | 4.30  | 5.37  | 14.47 | 12.90 | 12.89 |
| HTZ1    | 0.87 | 1.77  | 1.36  | 1.92 | 2.83  | 4.35  | 9.85  | 22.07 | 42.25 |
| YLR349W | 0.80 | 2.59  | 4.10  | 2.81 | 3.47  | 5.12  | 8.47  | 15.73 | 16.03 |
| STE50   | 0.84 | 3.74  | 2.35  | 3.75 | 2.95  | 7.85  | 9.81  | 11.66 | 15.86 |
| MAD1    | 0.88 | 1.21  | 1.26  | 1.44 | 2.18  | 5.69  | 8.87  | 21.32 | 24.96 |
| PHO80   | 1.15 | 0.89  | 1.13  | 1.44 | 1.51  | 8.90  | 8.65  | 19.05 | 25.57 |
| ALF1    | 0.75 | 0.14  | 0.27  | 1.63 | 4.96  | 18.95 | 6.65  | 8.64  | 6.84  |
| MAE1    | 0.85 | 6.65  | 7.20  | 4.34 | 5.01  | 3.03  | 6.59  | 8.00  | 23.63 |
| VID21   | 0.98 | 3.36  | 1.92  | 2.90 | 6.66  | 3.11  | 13.21 | 8.83  | 9.35  |
| DIC1    | 1.15 | 4.11  | 5.27  | 3.97 | 5.29  | 3.84  | 8.76  | 7.98  | 8.54  |
| VPS34   | 0.91 | 1.77  | 1.30  | 1.84 | 6.30  | 9.16  | 4.35  | 14.52 | 8.96  |
| MRPS8   | 0.64 | 1.32  | 1.06  | 1.00 | 3.35  | 7.64  | 8.97  | 15.85 | 14.21 |
| MRPL16  | 0.95 | 3.58  | 2.40  | 1.08 | 1.98  | 12.68 | 8.40  | 8.23  | 7.37  |
| GUF1    | 1.23 | 4.31  | 3.93  | 2.79 | 2.29  | 5.01  | 6.42  | 12.75 | 18.83 |
| YCR006C | 0.56 | 1.09  | 6.30  | 2.15 | 1.55  | 2.87  | 18.71 | 5.35  | 4.26  |
| PIN4    | 1.03 | 2.06  | 1.87  | 1.31 | 2.13  | 4.14  | 9.95  | 14.71 | 21.54 |
| BUB1    | 1.12 | 1.58  | 0.88  | 5.34 | 4.96  | 5.01  | 8.46  | 9.81  | 6.69  |
| TIM18   | 0.85 | 2.86  | 2.68  | 5.27 | 10.33 | 6.18  | 3.86  | 5.12  | 3.07  |
| FZO1    | 1.13 | 2.47  | 2.34  | 3.36 | 3.84  | 5.77  | 15.03 | 3.00  | 46.29 |
| BUB3    | 1.01 | 3.22  | 2.67  | 3.82 | 4.24  | 8.85  | 5.11  | 7.87  | 4.32  |
| UBP6    | 1.22 | 0.82  | 0.58  | 2.01 | 4.50  | 9.70  | 12.38 | 5.11  | 4.74  |
| YML090W | 1.21 | 2.20  | 2.99  | 2.18 | 2.54  | 5.57  | 9.46  | 10.06 | 16.73 |
| VAC7    | 0.92 | 0.47  | 0.46  | 2.16 | 2.22  | 19.60 | 5.50  | 4.70  | 11.83 |
| ATP12   | 1.00 | 3.96  | 3.58  | 3.01 | 3.42  | 6.26  | 4.98  | 9.43  | 9.68  |
| YLR296W | 1.20 | 4.09  | 3.26  | 2.29 | 1.93  | 4.13  | 7.63  | 10.38 | 3.60  |
| ADO1    | 1.77 | 1.69  | 0.26  | 1.38 | 1.16  | 17.90 | 6.18  | 4.44  | 8.32  |
| YCR079W | 0.88 | 3.68  | 3.22  | 3.26 | 5.15  | 8.74  | 6.03  | 3.82  | 2.80  |
| YHR100C | 1.10 | 2.76  | 2.19  | 2.23 | 2.40  | 5.61  | 8.13  | 10.25 | 40.40 |
| MDM30   | 1.26 | 1.92  | 2.50  | 3.09 | 6.99  | 4.78  | 6.50  | 7.58  | 10.06 |

|         |      |      |      |      |      |       |       |       |       |
|---------|------|------|------|------|------|-------|-------|-------|-------|
| MDM38   | 1.05 | 1.71 | 2.30 | 1.69 | 1.78 | 10.78 | 6.44  | 8.64  | 10.96 |
| RPS6B   | 1.66 | 1.31 | 0.99 | 1.20 | 1.03 | 15.76 | 6.65  | 5.33  | 6.75  |
| SIR3    | 0.80 | 1.89 | 1.00 | 1.21 | 2.62 | 4.97  | 14.43 | 6.90  | 23.29 |
| STE20   | 1.14 | 1.41 | 2.01 | 2.77 | 2.71 | 4.49  | 10.77 | 8.44  | 34.24 |
| MIS1    | 1.01 | 2.09 | 5.34 | 2.64 | 3.50 | 5.48  | 4.89  | 8.70  | 18.80 |
| NOT3    | 1.04 | 2.46 | 3.58 | 3.98 | 4.86 | 7.81  | 4.25  | 5.57  | 11.66 |
| YAF9    | 0.75 | 2.29 | 1.47 | 1.71 | 1.96 | 4.95  | 5.78  | 14.54 | 34.00 |
| YER087W | 1.10 | 2.57 | 1.49 | 1.17 | 1.70 | 3.85  | 6.19  | 14.57 | 16.68 |
| EMI1    | 1.02 | 3.96 | 3.34 | 3.19 | 3.86 | 3.25  | 6.59  | 7.25  | 8.09  |
| MEC3    | 1.00 | 3.00 | 2.95 | 2.01 | 2.24 | 5.54  | 7.08  | 8.57  | 7.98  |
| ARP6    | 0.97 | 1.11 | 1.07 | 1.65 | 2.27 | 3.32  | 7.61  | 14.26 | 30.10 |
| MRP51   | 1.25 | 2.92 | 1.05 | 1.12 | 2.18 | 3.46  | 7.99  | 12.16 | 6.02  |
| YKL030W | 1.01 | 4.05 | 5.14 | 5.67 | 5.83 | 2.72  | 3.99  | 3.67  | 4.58  |
| YGL196W | 1.17 | 0.89 | 1.14 | 0.93 | 1.21 | 4.25  | 6.04  | 16.44 | 13.87 |
| SOV1    | 1.16 | 2.87 | 1.91 | 1.37 | 3.03 | 5.43  | 8.45  | 7.73  | 6.38  |
| ATP11   | 0.97 | 4.35 | 3.19 | 1.69 | 3.00 | 2.96  | 5.38  | 10.33 | 9.05  |
| YLR290C | 0.99 | 4.24 | 2.77 | 2.51 | 1.84 | 5.16  | 8.51  | 5.79  | 5.87  |
| SIR3    | 0.91 | 2.10 | 1.15 | 1.58 | 2.06 | 8.19  | 9.21  | 5.73  | 17.42 |
| YLR257W | 1.02 | 1.51 | 1.51 | 1.26 | 1.37 | 6.12  | 9.53  | 8.53  | 20.99 |
| CSR2    | 1.16 | 1.66 | 1.00 | 1.06 | 1.09 | 6.16  | 5.30  | 13.40 | 5.29  |
| YOR205C | 1.17 | 2.78 | 2.04 | 1.65 | 2.21 | 3.78  | 6.94  | 10.23 | 12.84 |
| BUB2    | 0.95 | 0.88 | 1.51 | 1.16 | 1.44 | 2.97  | 4.77  | 16.94 | 10.23 |
| MDH2    | 0.94 | 3.52 | 5.43 | 4.62 | 3.08 | 2.68  | 4.29  | 6.01  | 7.58  |
| RPO41   | 1.06 | 3.94 | 2.92 | 2.39 | 2.98 | 3.20  | 3.34  | 10.67 | 4.85  |
| SLT2    | 0.86 | 1.80 | 1.14 | 0.91 | 1.66 | 8.25  | 8.43  | 7.26  | 3.37  |
| PAN3    | 1.12 | 0.76 | 1.28 | 1.81 | 2.25 | 3.68  | 5.76  | 13.64 | 30.16 |
| YLR280C | 1.03 | 4.19 | 2.85 | 2.31 | 1.79 | 4.13  | 7.63  | 6.30  | 4.61  |
| MSH2    | 0.69 | 2.85 | 3.33 | 1.98 | 2.70 | 4.63  | 8.10  | 5.82  | 14.11 |
| MSU1    | 0.78 | 2.60 | 1.46 | 1.23 | 1.67 | 3.83  | 5.05  | 12.97 | 7.56  |
| YPR099C | 1.09 | 3.33 | 1.20 | 1.33 | 2.43 | 4.72  | 7.43  | 8.00  | 6.97  |
| YLR312C | 1.06 | 4.21 | 3.36 | 2.13 | 2.77 | 4.65  | 4.88  | 6.33  | 10.10 |
| ATP5    | 0.87 | 5.45 | 4.30 | 1.92 | 2.73 | 4.53  | 4.49  | 5.03  | 5.36  |
| ADH3    | 0.83 | 1.43 | 2.00 | 1.93 | 3.59 | 5.59  | 4.66  | 8.97  | 6.86  |
| DEF1    | 2.91 | 0.29 | 0.62 | 0.48 | 1.27 | 4.91  | 3.19  | 15.19 | 4.52  |
| NAT1    | 1.02 | 1.32 | 0.72 | 1.61 | 2.67 | 3.64  | 6.38  | 11.35 | 19.54 |
| IMG1    | 1.33 | 3.19 | 2.71 | 1.59 | 1.95 | 5.15  | 6.67  | 6.11  | 10.21 |
| YJL022W | 0.86 | 2.43 | 0.78 | 0.58 | 1.04 | 2.19  | 4.29  | 16.34 | 11.76 |
| BCK1    | 0.82 | 1.82 | 1.14 | 1.26 | 1.56 | 7.48  | 7.30  | 6.81  | 2.72  |
| YLR283W | 1.16 | 4.23 | 3.14 | 2.71 | 2.09 | 4.39  | 3.70  | 6.73  | 3.61  |
| SDH2    | 1.05 | 2.40 | 3.32 | 1.29 | 1.66 | 3.89  | 4.65  | 9.73  | 4.26  |
| MRP1    | 0.99 | 3.82 | 1.19 | 1.30 | 2.18 | 3.56  | 7.08  | 7.66  | 11.28 |
| TCM10   | 0.82 | 3.13 | 2.75 | 1.95 | 3.25 | 2.43  | 7.28  | 6.02  | 3.47  |
| GGA1    | 1.01 | 1.18 | 0.79 | 0.78 | 1.07 | 4.04  | 12.41 | 6.23  | 21.82 |
| PET8    | 0.66 | 2.74 | 1.47 | 1.26 | 2.07 | 2.70  | 3.73  | 12.87 | 7.02  |

|           |      |      |      |      |      |       |       |       |       |
|-----------|------|------|------|------|------|-------|-------|-------|-------|
| YLR287C   | 1.14 | 3.02 | 2.50 | 1.44 | 2.25 | 3.55  | 3.36  | 10.22 | 2.83  |
| ATP10     | 0.86 | 1.84 | 2.77 | 1.53 | 1.86 | 2.73  | 6.30  | 9.52  | 5.64  |
| YLR279W   | 1.09 | 2.83 | 2.86 | 2.06 | 2.09 | 4.92  | 3.92  | 7.61  | 3.94  |
| MGM101    | 0.98 | 4.01 | 2.02 | 1.85 | 2.54 | 3.22  | 4.46  | 8.27  | 13.22 |
| ATP7      | 1.09 | 4.06 | 2.87 | 2.94 | 2.70 | 4.17  | 5.20  | 4.28  | 4.54  |
| SHP1      | 1.14 | 1.25 | 1.25 | 0.88 | 1.07 | 4.57  | 4.02  | 13.10 | 3.39  |
| SNF8      | 1.09 | 1.76 | 1.14 | 4.39 | 4.65 | 7.92  | 3.03  | 3.12  | 6.06  |
| SLS1      | 0.88 | 1.04 | 1.06 | 1.64 | 2.39 | 2.33  | 12.01 | 5.65  | 19.27 |
| DIA4      | 1.41 | 2.03 | 2.77 | 1.70 | 2.01 | 2.64  | 5.13  | 9.24  | 4.34  |
| MSE1      | 1.26 | 1.36 | 2.46 | 1.02 | 1.55 | 2.70  | 9.16  | 7.17  | 5.09  |
| MDM12     | 1.22 | 1.50 | 1.19 | 0.72 | 0.93 | 3.05  | 3.31  | 14.76 | 5.68  |
| EXG1      | 1.15 | 2.26 | 2.47 | 1.46 | 1.80 | 4.02  | 3.74  | 9.67  | 8.89  |
| POS5      | 1.05 | 2.76 | 4.96 | 1.86 | 0.64 | 6.85  | 3.72  | 4.72  | 13.02 |
| PDC1      | 0.92 | 1.27 | 2.17 | 1.79 | 1.94 | 8.18  | 6.60  | 3.67  | 5.37  |
| YLR358C   | 1.13 | 2.11 | 2.83 | 1.78 | 3.77 | 5.24  | 3.57  | 5.95  | 8.60  |
| ECM38     | 1.22 | 3.51 | 3.21 | 1.67 | 1.69 | 4.01  | 7.18  | 3.86  | 2.44  |
| MTF1      | 0.75 | 1.72 | 1.53 | 1.51 | 2.26 | 4.45  | 6.78  | 7.25  | 18.05 |
| GRX5      | 0.97 | 1.28 | 0.98 | 0.92 | 1.61 | 7.12  | 5.59  | 7.67  | 5.21  |
| YPL005W   | 0.67 | 3.79 | 2.52 | 1.72 | 2.24 | 2.59  | 4.66  | 7.96  | 3.75  |
| LPE10     | 1.16 | 1.89 | 1.87 | 1.95 | 2.17 | 4.26  | 4.55  | 8.27  | 5.99  |
| TOM70     | 0.78 | 1.55 | 1.83 | 1.85 | 2.93 | 4.42  | 5.93  | 6.83  | 15.25 |
| WSC3      | 0.99 | 2.25 | 4.13 | 1.95 | 1.59 | 3.83  | 7.63  | 3.73  | 7.77  |
| URE2      | 0.76 | 1.16 | 3.94 | 2.35 | 3.47 | 3.29  | 4.48  | 6.58  | 8.59  |
| MRP2      | 0.80 | 2.36 | 1.77 | 0.66 | 2.18 | 2.22  | 6.94  | 9.03  | 8.37  |
| YJR120W   | 1.05 | 3.75 | 3.24 | 2.54 | 2.26 | 3.79  | 3.74  | 5.55  | 4.24  |
| YHR039C-B | 2.01 | 0.21 | 0.44 | 1.45 | 1.03 | 11.14 | 2.71  | 6.91  | 6.81  |
| MRPL32    | 0.78 | 4.06 | 1.40 | 1.10 | 2.49 | 2.39  | 4.36  | 9.28  | 8.81  |
| MRPL22    | 1.01 | 2.67 | 0.99 | 1.03 | 1.63 | 2.06  | 4.35  | 12.09 | 5.31  |
| YDL033C   | 1.15 | 1.58 | 1.55 | 1.78 | 2.21 | 3.64  | 5.62  | 8.19  | 9.83  |
| MSM1      | 1.30 | 3.64 | 1.75 | 1.85 | 3.03 | 3.57  | 4.57  | 5.86  | 9.66  |
| YER077C   | 0.97 | 1.99 | 1.14 | 0.98 | 2.34 | 1.88  | 5.26  | 10.97 | 13.66 |
| BIK1      | 0.90 | 1.54 | 1.61 | 1.62 | 2.03 | 4.72  | 7.68  | 5.42  | 18.38 |
| BFA1      | 1.22 | 0.88 | 1.12 | 1.11 | 1.23 | 3.17  | 5.58  | 11.22 | 11.48 |
| VMA5      | 0.89 | 0.12 | 0.31 | 0.89 | 1.22 | 13.27 | 5.31  | 3.37  | 10.25 |
| MRPL33    | 0.94 | 3.55 | 1.43 | 1.26 | 1.84 | 2.92  | 8.43  | 4.87  | 3.33  |
| NPR2      | 0.88 | 1.18 | 1.79 | 2.07 | 3.24 | 2.80  | 6.11  | 7.13  | 11.25 |
| YCL003W   | 1.05 | 1.63 | 1.41 | 1.24 | 1.61 | 3.81  | 2.61  | 11.77 | 7.32  |
| CDA2      | 1.22 | 3.32 | 3.94 | 1.88 | 1.75 | 3.91  | 3.71  | 5.30  | 10.27 |
| NKP2      | 0.97 | 3.42 | 2.84 | 1.73 | 2.17 | 5.98  | 4.12  | 3.76  | 2.31  |
| YJL169W   | 0.95 | 3.08 | 1.28 | 1.50 | 1.90 | 2.26  | 7.75  | 6.24  | 5.92  |
| YPL107W   | 1.13 | 3.19 | 2.69 | 2.13 | 1.45 | 3.70  | 4.60  | 6.05  | 4.31  |
| ATP11     | 0.67 | 2.89 | 1.76 | 1.31 | 1.67 | 2.39  | 4.33  | 9.91  | 8.90  |
| YPL105C   | 0.90 | 1.36 | 1.07 | 1.57 | 1.45 | 3.21  | 3.88  | 11.42 | 6.63  |
| KAP123    | 1.01 | 1.60 | 4.65 | 1.67 | 3.80 | 2.55  | 5.09  | 4.47  | 4.94  |



|           |      |      |      |      |      |      |      |       |       |
|-----------|------|------|------|------|------|------|------|-------|-------|
| RSM22     | 1.12 | 3.07 | 2.31 | 1.62 | 1.85 | 3.75 | 5.11 | 5.99  | 7.21  |
| YDL025C   | 1.02 | 1.78 | 3.01 | 3.52 | 5.58 | 3.34 | 3.04 | 3.52  | 6.48  |
| AEP2      | 0.94 | 3.19 | 2.54 | 1.99 | 2.00 | 3.25 | 4.51 | 6.32  | 9.48  |
| HCM1      | 0.79 | 2.19 | 1.41 | 1.37 | 1.87 | 3.37 | 4.37 | 9.35  | 25.53 |
| MRPL4     | 1.18 | 2.56 | 2.14 | 1.10 | 2.93 | 3.77 | 5.16 | 5.80  | 6.43  |
| PET56     | 0.72 | 2.06 | 1.18 | 1.11 | 1.38 | 2.38 | 5.07 | 10.72 | 11.61 |
| PET130    | 1.20 | 2.95 | 1.80 | 1.06 | 1.81 | 3.59 | 2.66 | 9.45  | 9.48  |
| MRPL9     | 1.03 | 2.21 | 1.12 | 1.18 | 1.61 | 2.43 | 3.50 | 11.22 | 7.12  |
| SET2      | 1.20 | 2.21 | 1.76 | 1.36 | 1.65 | 3.13 | 5.59 | 7.37  | 20.93 |
| NPR1      | 1.30 | 3.95 | 3.27 | 2.01 | 2.65 | 2.53 | 3.37 | 5.17  | 8.09  |
| MAD3      | 2.14 | 1.96 | 2.55 | 2.06 | 2.99 | 1.42 | 4.39 | 6.73  | 10.81 |
| MSK1      | 0.90 | 2.12 | 0.99 | 0.51 | 1.59 | 2.32 | 7.35 | 8.42  | 8.32  |
| EMI1      | 0.74 | 2.64 | 2.76 | 2.86 | 2.78 | 3.32 | 3.25 | 5.75  | 3.34  |
| YIL064W   | 1.09 | 2.91 | 2.73 | 2.20 | 2.41 | 3.85 | 3.27 | 5.59  | 2.81  |
| RSM19     | 0.65 | 3.03 | 0.79 | 1.63 | 3.26 | 2.29 | 3.02 | 9.34  | 3.74  |
| ISM1      | 1.22 | 2.44 | 2.25 | 1.28 | 1.35 | 2.50 | 4.27 | 8.69  | 6.09  |
| GLC8      | 0.96 | 1.63 | 1.40 | 2.10 | 2.31 | 4.56 | 5.60 | 5.39  | 8.92  |
| CDA1      | 1.06 | 1.92 | 2.29 | 1.41 | 1.25 | 3.11 | 5.91 | 6.99  | 7.78  |
| YLR281C   | 0.86 | 3.80 | 2.98 | 1.84 | 1.60 | 3.65 | 5.94 | 3.20  | 14.11 |
| YGR150C   | 0.75 | 2.17 | 1.71 | 0.95 | 2.02 | 4.13 | 5.38 | 6.72  | 4.22  |
| MSC1      | 0.80 | 3.52 | 2.75 | 1.72 | 1.75 | 4.65 | 2.86 | 5.73  | 22.98 |
| AOR1      | 0.67 | 1.11 | 0.86 | 1.11 | 2.26 | 2.45 | 2.76 | 12.49 | 20.07 |
| RMD12     | 1.25 | 1.64 | 1.18 | 1.12 | 2.43 | 5.14 | 5.65 | 5.30  | 5.66  |
| MCK1      | 1.02 | 1.72 | 1.79 | 1.24 | 1.79 | 4.57 | 5.08 | 6.15  | 4.06  |
| PDR15     | 1.00 | 0.62 | 0.61 | 0.61 | 0.68 | 7.79 | 6.73 | 5.29  | 6.62  |
| YDR065W   | 1.28 | 2.54 | 2.02 | 1.21 | 1.78 | 3.75 | 3.43 | 7.33  | 8.27  |
| MRPL49    | 0.90 | 1.87 | 0.79 | 0.75 | 1.58 | 1.96 | 6.14 | 9.10  | 10.35 |
| YPL098C   | 1.15 | 2.65 | 1.62 | 1.57 | 2.09 | 3.38 | 4.63 | 6.01  | 3.97  |
| YNL213C   | 0.86 | 2.93 | 1.76 | 1.03 | 1.47 | 2.90 | 5.95 | 6.19  | 13.59 |
| YMR293C   | 0.74 | 2.46 | 1.04 | 1.20 | 1.60 | 4.10 | 3.59 | 8.18  | 10.58 |
| SPH1      | 1.34 | 3.49 | 3.42 | 2.29 | 1.85 | 3.10 | 3.98 | 3.27  | 6.23  |
| SWC1      | 1.15 | 1.41 | 1.26 | 3.52 | 2.40 | 3.76 | 4.55 | 4.58  | 12.18 |
| MRH4      | 0.94 | 1.73 | 1.29 | 1.33 | 1.61 | 2.66 | 4.93 | 8.13  | 79.01 |
| RSC1      | 0.81 | 1.68 | 3.24 | 3.28 | 4.85 | 2.76 | 2.80 | 3.13  | 4.13  |
| ATP4      | 0.97 | 3.08 | 2.79 | 2.26 | 1.59 | 3.76 | 2.35 | 5.74  | 4.80  |
| UBC12     | 1.09 | 2.79 | 3.45 | 1.59 | 1.99 | 2.81 | 4.11 | 4.66  | 8.71  |
| YPL183W-A | 1.15 | 2.87 | 1.93 | 1.51 | 1.82 | 2.79 | 4.29 | 6.08  | 3.13  |
| YDR532C   | 0.98 | 1.88 | 0.60 | 0.88 | 2.00 | 3.65 | 7.78 | 4.65  | 3.31  |
| MRPL23    | 1.19 | 3.37 | 1.44 | 1.27 | 1.09 | 3.31 | 4.53 | 6.19  | 13.00 |
| HOS2      | 0.96 | 1.84 | 0.93 | 1.61 | 2.39 | 2.23 | 2.60 | 9.80  | 5.36  |
| YKL169C   | 0.73 | 2.43 | 1.37 | 1.38 | 1.17 | 5.79 | 3.76 | 5.72  | 5.57  |
| MDM10     | 0.82 | 2.44 | 0.74 | 1.11 | 1.00 | 1.08 | 2.37 | 12.72 | 10.66 |
| STP1      | 0.87 | 2.63 | 4.45 | 2.64 | 2.47 | 2.28 | 3.17 | 3.77  | 2.40  |
| YLL033W   | 0.94 | 1.80 | 1.52 | 1.31 | 2.06 | 3.72 | 3.24 | 7.68  | 23.79 |

|           |      |      |      |      |      |      |      |       |       |
|-----------|------|------|------|------|------|------|------|-------|-------|
| MRPS9     | 0.81 | 1.25 | 0.97 | 1.01 | 2.21 | 3.31 | 4.91 | 7.64  | 4.91  |
| YPR116W   | 0.48 | 1.37 | 1.14 | 1.14 | 1.75 | 2.69 | 3.41 | 10.13 | 6.13  |
| PSR1      | 0.72 | 2.58 | 1.79 | 1.23 | 1.51 | 3.88 | 2.46 | 7.93  | 8.94  |
| VIP1      | 1.18 | 2.36 | 2.65 | 1.71 | 1.93 | 4.33 | 4.63 | 3.29  | 4.85  |
| MIP1      | 1.26 | 1.65 | 1.76 | 1.38 | 3.33 | 3.73 | 5.64 | 3.33  | 9.34  |
| ARG80     | 0.65 | 2.94 | 1.10 | 1.05 | 1.83 | 4.40 | 3.61 | 6.49  | 7.85  |
| MRPS16    | 1.05 | 2.37 | 1.37 | 1.16 | 2.09 | 2.64 | 6.23 | 5.12  | 5.44  |
| RSM18     | 0.73 | 3.92 | 1.25 | 0.95 | 2.06 | 4.63 | 2.87 | 5.62  | 11.55 |
| MRPS17    | 1.01 | 2.49 | 1.93 | 1.10 | 1.30 | 3.07 | 4.12 | 6.97  | 4.93  |
| YOR200W   | 0.85 | 3.33 | 1.57 | 1.15 | 1.87 | 2.23 | 4.39 | 6.53  | 10.81 |
| VPS36     | 1.21 | 1.07 | 1.41 | 3.27 | 4.71 | 3.72 | 3.75 | 2.66  | 5.68  |
| GCV1      | 1.07 | 1.55 | 1.77 | 1.81 | 2.51 | 2.92 | 3.66 | 6.45  | 3.47  |
| VMA8      | 0.90 | 0.16 | 0.26 | 0.84 | 1.61 | 3.15 | 2.11 | 12.65 | 7.30  |
| IMG2      | 0.69 | 0.95 | 0.80 | 0.72 | 1.45 | 5.35 | 6.20 | 5.45  | 5.47  |
| MRF1      | 0.70 | 2.37 | 1.72 | 1.34 | 1.75 | 2.50 | 2.96 | 8.21  | 7.64  |
| STP22     | 1.37 | 1.43 | 0.60 | 2.76 | 3.73 | 4.66 | 4.30 | 2.67  | 4.45  |
| MRPL19    | 0.70 | 1.76 | 3.40 | 1.01 | 1.40 | 5.02 | 3.11 | 5.11  | 16.12 |
| QRI5      | 1.16 | 2.65 | 1.93 | 1.57 | 3.11 | 3.79 | 3.48 | 3.77  | 12.59 |
| MSY1      | 0.84 | 1.43 | 0.84 | 0.98 | 1.48 | 2.86 | 7.17 | 5.87  | 12.53 |
| VPS25     | 1.26 | 1.37 | 1.27 | 2.70 | 3.80 | 4.39 | 3.78 | 2.88  | 3.70  |
| IFM1      | 0.65 | 2.95 | 1.52 | 1.31 | 2.82 | 3.80 | 2.83 | 5.54  | 10.14 |
| YGR102C   | 0.88 | 2.66 | 1.68 | 0.89 | 1.50 | 3.78 | 3.65 | 6.28  | 10.85 |
| MRP13     | 1.16 | 2.40 | 1.79 | 1.53 | 2.15 | 3.16 | 5.74 | 3.37  | 10.38 |
| ARD1      | 1.26 | 1.72 | 0.69 | 1.18 | 1.89 | 2.38 | 3.16 | 8.93  | 18.94 |
| BUD6      | 0.85 | 2.56 | 2.40 | 1.70 | 1.45 | 3.71 | 4.04 | 4.40  | 1.67  |
| YDR479C   | 1.00 | 1.92 | 1.89 | 2.67 | 2.89 | 2.99 | 3.57 | 4.14  | 5.53  |
| MTG1      | 0.90 | 2.02 | 0.94 | 0.94 | 1.29 | 2.90 | 4.15 | 7.91  | 14.68 |
| YDR248C   | 0.87 | 1.05 | 0.99 | 0.95 | 0.69 | 5.45 | 6.72 | 4.19  | 94.77 |
| PDE2      | 0.95 | 1.89 | 1.37 | 0.85 | 1.34 | 7.37 | 4.28 | 2.84  | 6.37  |
| YIL077C   | 1.20 | 3.16 | 2.70 | 1.96 | 2.22 | 3.40 | 2.52 | 3.57  | 5.02  |
| PUS2      | 0.96 | 1.69 | 2.36 | 1.87 | 1.90 | 2.63 | 2.78 | 6.48  | 4.92  |
| AEP1      | 0.85 | 2.46 | 1.77 | 1.24 | 1.76 | 2.20 | 4.65 | 5.70  | 5.94  |
| YLR282C   | 0.92 | 3.02 | 2.50 | 1.83 | 1.99 | 4.15 | 3.66 | 2.51  | 3.86  |
| CAJ1      | 1.00 | 2.28 | 1.53 | 1.72 | 2.24 | 2.81 | 3.85 | 5.15  | 5.60  |
| YDL032W   | 1.12 | 1.96 | 1.70 | 1.16 | 1.50 | 2.73 | 4.29 | 6.10  | 6.40  |
| MRP20     | 1.03 | 3.01 | 0.87 | 1.44 | 1.74 | 3.39 | 3.04 | 6.03  | 7.42  |
| YPL099C   | 1.12 | 1.88 | 2.26 | 1.65 | 2.12 | 2.48 | 2.81 | 6.14  | 4.12  |
| MON2      | 0.86 | 1.86 | 2.05 | 1.26 | 1.24 | 2.58 | 4.33 | 6.24  | 19.23 |
| PET123    | 1.29 | 1.93 | 1.77 | 1.37 | 1.75 | 3.53 | 3.23 | 5.50  | 5.04  |
| YML100W-A | 0.76 | 2.47 | 1.40 | 1.20 | 1.85 | 4.06 | 4.33 | 3.99  | 19.45 |
| THR4      | 1.15 | 1.11 | 2.49 | 3.22 | 2.38 | 3.13 | 3.55 | 3.01  | 5.01  |
| SUV3      | 0.97 | 2.22 | 1.35 | 1.00 | 1.19 | 2.89 | 4.75 | 5.63  | 8.92  |
| PIF1      | 1.05 | 2.19 | 2.81 | 1.94 | 3.45 | 1.77 | 2.21 | 4.55  | 8.27  |
| YBR004C   | 0.96 | 1.71 | 0.89 | 0.84 | 1.45 | 1.71 | 2.43 | 9.94  | 15.60 |

|         |      |      |      |      |      |      |      |      |       |
|---------|------|------|------|------|------|------|------|------|-------|
| YDR114C | 1.09 | 1.27 | 1.93 | 1.14 | 2.14 | 3.00 | 3.67 | 5.66 | 10.17 |
| YLR091W | 1.04 | 2.58 | 1.47 | 1.12 | 1.40 | 2.74 | 3.00 | 6.49 | 5.54  |
| SET3    | 0.78 | 1.81 | 0.63 | 1.47 | 1.91 | 3.35 | 2.98 | 6.86 | 11.36 |
| SNQ2    | 1.08 | 0.78 | 0.56 | 0.88 | 1.02 | 2.61 | 5.13 | 7.67 | 25.47 |
| YLR311C | 0.78 | 3.21 | 2.76 | 1.55 | 1.68 | 1.90 | 3.17 | 4.68 | 7.73  |
| MRPL8   | 1.22 | 1.96 | 1.22 | 1.17 | 1.99 | 2.12 | 5.02 | 4.98 | 4.45  |
| EAF3    | 1.01 | 2.31 | 1.61 | 1.33 | 1.41 | 1.81 | 4.54 | 5.47 | 12.94 |
| ATP2    | 1.00 | 2.54 | 3.50 | 1.94 | 2.02 | 3.10 | 2.00 | 3.38 | 7.85  |
| NNT1    | 1.27 | 2.66 | 2.60 | 1.39 | 1.66 | 3.61 | 2.58 | 3.68 | 5.31  |
| AQR1    | 0.91 | 1.74 | 2.22 | 2.89 | 2.54 | 1.82 | 3.15 | 4.14 | 3.55  |
| YLF2    | 0.70 | 3.15 | 0.95 | 0.81 | 2.03 | 3.02 | 4.06 | 4.62 | 3.33  |
| RML2    | 1.27 | 2.24 | 1.30 | 0.98 | 1.68 | 3.33 | 4.93 | 3.60 | 7.78  |
| YJL045W | 0.91 | 1.46 | 2.94 | 1.76 | 2.04 | 2.74 | 3.38 | 4.08 | 3.28  |
| LTE1    | 0.96 | 1.29 | 3.01 | 3.56 | 3.46 | 1.67 | 2.00 | 3.31 | 6.52  |
| YPL184C | 1.09 | 1.98 | 1.08 | 1.05 | 1.26 | 2.14 | 4.97 | 5.69 | 9.30  |
| PRP12   | 1.52 | 1.03 | 0.04 | 0.71 | 1.00 | 5.93 | 5.14 | 3.82 | 3.93  |
| REF2    | 0.97 | 0.82 | 1.55 | 1.94 | 1.80 | 5.45 | 3.85 | 2.81 | 7.50  |
| ECI1    | 0.90 | 2.32 | 2.58 | 1.69 | 1.47 | 3.50 | 2.89 | 3.84 | 13.29 |
| MGM1    | 1.13 | 1.36 | 1.82 | 1.43 | 2.56 | 4.41 | 2.61 | 3.81 | 5.30  |
| OAR1    | 1.29 | 1.03 | 1.36 | 1.55 | 1.79 | 3.12 | 4.29 | 4.69 | 8.11  |
| MRPL24  | 1.31 | 1.65 | 1.65 | 1.35 | 2.05 | 2.46 | 3.60 | 5.00 | 7.93  |
| PET112  | 0.63 | 1.30 | 0.82 | 1.04 | 1.08 | 3.95 | 4.88 | 5.33 | 11.12 |
| YBR219C | 0.74 | 2.52 | 2.14 | 2.23 | 2.97 | 2.19 | 2.61 | 3.60 | 6.87  |
| MAK10   | 1.07 | 1.85 | 1.76 | 1.25 | 1.18 | 1.82 | 2.83 | 7.23 | 7.21  |
| YBR281C | 1.19 | 3.21 | 1.68 | 1.65 | 1.92 | 3.64 | 1.71 | 3.94 | 3.78  |
| MDM20   | 0.79 | 1.31 | 0.80 | 1.43 | 2.06 | 2.08 | 4.22 | 6.04 | 16.78 |
| BUD26   | 1.13 | 0.77 | 0.51 | 0.83 | 1.06 | 1.93 | 4.42 | 8.06 | 34.76 |
| AHP1    | 0.82 | 0.89 | 0.91 | 0.92 | 1.07 | 3.43 | 2.89 | 7.69 | 39.28 |
| SUM1    | 1.32 | 1.27 | 2.34 | 2.18 | 3.06 | 2.04 | 2.49 | 3.78 | 4.98  |
| MSN1    | 0.97 | 2.21 | 1.65 | 1.69 | 2.46 | 2.37 | 3.06 | 4.00 | 3.05  |
| RSM27   | 1.07 | 2.44 | 1.29 | 1.10 | 2.36 | 2.61 | 4.54 | 2.82 | 2.36  |
| MRPL38  | 1.11 | 2.33 | 1.53 | 0.87 | 1.11 | 2.36 | 3.46 | 5.35 | 3.09  |
| MRS4    | 1.17 | 1.65 | 1.47 | 1.18 | 1.32 | 2.59 | 3.15 | 5.53 | 1.75  |
| MST1    | 0.74 | 2.98 | 1.56 | 1.00 | 1.39 | 1.74 | 4.27 | 4.31 | 4.45  |
| YCR050C | 1.17 | 2.10 | 1.97 | 1.01 | 1.24 | 4.57 | 1.75 | 4.17 | 7.96  |
| SPC72   | 1.13 | 1.62 | 0.78 | 1.85 | 2.69 | 3.33 | 2.71 | 3.73 | 12.37 |
| YJL149W | 0.85 | 1.82 | 2.14 | 1.84 | 1.84 | 2.12 | 2.94 | 4.21 | 4.99  |
| BUD25   | 1.04 | 1.55 | 1.28 | 1.36 | 1.66 | 2.51 | 3.24 | 4.98 | 0.86  |
| MIG1    | 0.56 | 1.56 | 2.60 | 1.03 | 1.62 | 3.97 | 2.93 | 3.35 | 6.36  |
| MRPL35  | 0.92 | 2.27 | 0.98 | 0.98 | 2.15 | 2.82 | 2.84 | 4.65 | 3.78  |
| YLR252W | 1.16 | 1.67 | 1.43 | 1.48 | 2.23 | 2.48 | 4.29 | 2.81 | 2.90  |
| MHR1    | 0.77 | 3.30 | 0.98 | 1.04 | 1.30 | 2.59 | 2.27 | 5.28 | 4.91  |
| NGR1    | 0.92 | 2.03 | 1.78 | 1.45 | 1.19 | 2.23 | 2.74 | 5.17 | 5.79  |
| PAK1    | 0.92 | 2.14 | 1.69 | 1.10 | 1.12 | 3.21 | 3.65 | 3.63 | 3.82  |

|           |      |      |      |      |      |      |      |      |       |
|-----------|------|------|------|------|------|------|------|------|-------|
| MRPL40    | 0.86 | 3.20 | 1.36 | 1.32 | 1.60 | 2.77 | 3.00 | 3.31 | 6.81  |
| MTF2      | 0.78 | 2.60 | 0.66 | 0.72 | 1.43 | 2.23 | 4.52 | 4.49 | 11.52 |
| NAM2      | 0.87 | 2.04 | 1.46 | 1.16 | 1.38 | 2.04 | 4.44 | 4.02 | 7.76  |
| HOG1      | 1.00 | 1.40 | 1.92 | 2.03 | 2.58 | 2.11 | 3.69 | 2.65 | 3.34  |
| KTR2      | 0.64 | 1.59 | 0.87 | 0.76 | 1.12 | 2.34 | 3.75 | 6.24 | 14.79 |
| SIR1      | 0.84 | 1.42 | 0.75 | 1.18 | 1.24 | 3.01 | 5.06 | 3.75 | 3.38  |
| RRN10     | 1.23 | 1.44 | 1.23 | 0.79 | 1.28 | 2.85 | 2.13 | 6.27 | 10.21 |
| HEM14     | 0.76 | 2.08 | 1.55 | 1.38 | 2.16 | 3.35 | 2.63 | 3.28 | 5.25  |
| RIM1      | 1.23 | 1.64 | 1.43 | 1.23 | 1.87 | 1.82 | 3.81 | 4.16 | 6.87  |
| SMY2      | 0.89 | 1.34 | 1.44 | 1.54 | 2.08 | 1.68 | 2.40 | 5.77 | 4.70  |
| YNL157W   | 0.74 | 0.78 | 0.65 | 0.66 | 0.82 | 6.45 | 3.00 | 3.99 | 7.19  |
| YOR199W   | 0.84 | 2.18 | 1.72 | 1.35 | 1.16 | 1.91 | 3.08 | 4.79 | 7.79  |
| SEC22     | 0.75 | 1.68 | 1.34 | 1.47 | 1.20 | 3.68 | 3.22 | 3.67 | 5.99  |
| RTN1      | 1.01 | 0.77 | 0.91 | 0.69 | 0.96 | 1.74 | 1.10 | 9.70 | 58.30 |
| SAM1      | 0.90 | 1.37 | 1.96 | 1.72 | 2.10 | 1.30 | 2.76 | 4.70 | 5.34  |
| VPS15     | 1.07 | 0.59 | 0.88 | 1.53 | 1.19 | 2.12 | 4.62 | 4.79 | 10.08 |
| YPL114W   | 1.08 | 2.35 | 2.13 | 1.47 | 1.57 | 2.69 | 2.89 | 2.57 | 2.27  |
| YLR235C   | 1.28 | 2.05 | 1.73 | 1.43 | 1.64 | 2.30 | 2.76 | 3.54 | 7.81  |
| SHM1      | 0.87 | 0.89 | 1.66 | 1.21 | 1.17 | 4.23 | 4.02 | 2.62 | 4.68  |
| NCL1      | 0.92 | 1.30 | 0.75 | 1.00 | 1.28 | 3.00 | 4.25 | 4.11 | 4.89  |
| HOS3      | 0.90 | 1.36 | 1.24 | 1.21 | 1.19 | 2.09 | 3.31 | 5.31 | 3.58  |
| YPL183W-A | 1.09 | 2.00 | 1.93 | 1.50 | 1.13 | 2.67 | 3.43 | 2.82 | 1.98  |
| LRS4      | 1.20 | 1.16 | 0.92 | 0.92 | 0.79 | 2.06 | 2.67 | 6.81 | 6.15  |
| YGL072C   | 1.02 | 3.53 | 0.86 | 0.75 | 1.58 | 1.74 | 3.41 | 3.54 | 12.68 |
| YNL296W   | 0.99 | 1.31 | 1.46 | 1.34 | 1.87 | 2.98 | 3.13 | 3.34 | 9.14  |
| ICY2      | 0.94 | 0.95 | 1.05 | 1.06 | 1.57 | 3.22 | 3.56 | 4.07 | 3.58  |
| TSA1      | 1.07 | 0.90 | 0.74 | 0.90 | 1.10 | 4.82 | 3.05 | 3.75 | 5.09  |
| MKK1      | 0.80 | 1.14 | 0.97 | 1.14 | 1.05 | 2.61 | 2.12 | 6.42 | 15.96 |
| SHE9      | 1.32 | 0.91 | 1.13 | 0.76 | 1.03 | 3.10 | 2.55 | 5.41 | 3.68  |
| YER093C-A | 0.83 | 2.49 | 1.99 | 1.95 | 1.72 | 2.13 | 2.18 | 2.88 | 3.83  |
| TTR1      | 0.83 | 1.34 | 1.45 | 1.60 | 1.93 | 2.82 | 2.37 | 3.77 | 7.67  |
| HMI1      | 0.88 | 1.47 | 1.85 | 1.31 | 1.27 | 2.30 | 3.27 | 3.71 | 11.60 |
| YLL029W   | 1.06 | 1.50 | 1.73 | 1.32 | 1.47 | 2.71 | 3.16 | 2.99 | 5.01  |
| YPL113C   | 1.07 | 2.12 | 1.51 | 1.24 | 1.47 | 1.98 | 3.25 | 3.30 | 2.79  |
| THI4      | 0.79 | 2.17 | 1.60 | 1.22 | 1.66 | 2.35 | 2.31 | 3.81 | 11.50 |
| SCD6      | 0.96 | 1.25 | 1.73 | 1.16 | 1.17 | 2.32 | 2.50 | 4.74 | 4.38  |
| YCR087C-A | 0.84 | 2.01 | 1.25 | 1.15 | 1.17 | 1.30 | 2.63 | 5.46 | 2.04  |
| MRPL11    | 0.94 | 2.51 | 1.20 | 1.41 | 0.78 | 2.52 | 2.96 | 3.48 | 3.73  |
| DEM1      | 1.15 | 2.37 | 1.34 | 1.51 | 1.46 | 1.87 | 1.98 | 4.04 | 5.55  |
| YMR075C-A | 1.11 | 2.06 | 2.07 | 1.09 | 1.20 | 1.74 | 2.51 | 3.93 | 9.65  |
| SWF5      | 1.44 | 1.19 | 1.32 | 1.18 | 1.24 | 3.11 | 2.50 | 3.69 | 6.28  |
| SMF1      | 0.73 | 1.54 | 1.61 | 1.12 | 1.30 | 2.11 | 3.42 | 3.81 | 3.43  |
| CCE1      | 1.28 | 0.93 | 1.64 | 1.41 | 1.74 | 2.29 | 3.52 | 2.80 | 2.36  |
| ETR1      | 0.78 | 1.39 | 1.17 | 1.12 | 1.25 | 2.61 | 3.30 | 3.86 | 3.00  |

|         |      |      |      |      |      |      |      |      |       |
|---------|------|------|------|------|------|------|------|------|-------|
| ICL1    | 1.23 | 2.84 | 0.75 | 1.33 | 1.85 | 2.13 | 1.82 | 3.53 | 5.50  |
| SSU1    | 1.01 | 1.24 | 1.42 | 1.32 | 1.71 | 1.82 | 2.69 | 4.27 | 4.45  |
| MRPL25  | 0.83 | 1.78 | 1.29 | 0.45 | 1.56 | 2.32 | 2.15 | 5.08 | 6.44  |
| SAS2    | 0.93 | 0.71 | 0.74 | 0.71 | 0.79 | 1.78 | 3.43 | 6.35 | 27.69 |
| MRPL37  | 1.02 | 2.49 | 0.89 | 0.85 | 1.80 | 2.24 | 3.38 | 2.74 | 5.01  |
| PMP3    | 1.07 | 0.72 | 0.52 | 0.99 | 0.75 | 2.71 | 2.08 | 6.56 | 8.20  |
| YPR050C | 0.89 | 1.45 | 1.33 | 1.14 | 1.11 | 1.71 | 2.62 | 5.07 | 7.42  |
| SEF1    | 1.13 | 2.05 | 1.27 | 1.16 | 1.21 | 1.49 | 2.14 | 4.88 | 4.39  |
| MAK31   | 0.66 | 1.90 | 1.23 | 1.45 | 1.10 | 1.38 | 2.80 | 4.77 | 4.16  |
| YLR199C | 0.67 | 0.64 | 0.60 | 0.74 | 0.97 | 2.44 | 5.91 | 3.28 | 14.74 |
| MSR1    | 1.09 | 1.93 | 0.93 | 0.98 | 1.57 | 1.21 | 3.35 | 4.15 | 17.47 |
| COX6    | 1.00 | 2.72 | 1.79 | 1.26 | 1.09 | 1.65 | 3.15 | 2.52 | 2.58  |
| HHO1    | 0.78 | 1.71 | 1.29 | 1.39 | 1.80 | 1.41 | 2.97 | 3.78 | 4.14  |
| COX19   | 1.12 | 1.49 | 1.30 | 1.33 | 1.32 | 1.99 | 3.02 | 3.40 | 4.34  |
| BUR2    | 0.50 | 0.53 | 0.77 | 1.83 | 0.98 | 2.76 | 2.69 | 4.70 | 2.66  |
| COX10   | 1.05 | 1.64 | 1.10 | 1.24 | 1.47 | 2.40 | 3.03 | 2.65 | 2.48  |
| CSM1    | 1.03 | 1.00 | 1.04 | 1.09 | 1.30 | 2.84 | 3.51 | 2.67 | 5.18  |
| YGK3    | 0.78 | 1.52 | 1.74 | 1.48 | 1.74 | 1.58 | 2.98 | 2.61 | 2.37  |
| YPL017C | 1.29 | 1.39 | 1.00 | 0.92 | 1.25 | 1.22 | 3.37 | 3.95 | 9.92  |
| IES6    | 1.01 | 1.41 | 0.97 | 1.23 | 1.20 | 1.97 | 3.78 | 2.77 | 6.88  |
| RBL2    | 0.90 | 0.73 | 1.10 | 1.51 | 2.04 | 1.94 | 2.65 | 3.46 | 33.40 |
| YOR052C | 1.01 | 1.44 | 0.95 | 1.11 | 1.34 | 1.85 | 2.37 | 4.07 | 4.55  |
| DCC1    | 0.95 | 1.48 | 1.34 | 0.92 | 1.28 | 1.78 | 3.06 | 3.34 | 16.03 |
| PET100  | 1.24 | 1.61 | 1.45 | 1.55 | 1.25 | 2.37 | 2.42 | 2.17 | 3.65  |
| ATE1    | 1.17 | 1.26 | 1.24 | 1.78 | 1.69 | 1.71 | 2.34 | 2.82 | 4.13  |
| STM1    | 1.09 | 1.34 | 0.66 | 1.17 | 1.38 | 2.58 | 2.48 | 3.27 | 11.56 |
| MRP49   | 1.05 | 1.70 | 1.52 | 1.45 | 1.25 | 1.63 | 2.40 | 2.91 | 4.33  |
| TRP1    | 1.10 | 1.60 | 0.83 | 1.69 | 1.47 | 1.75 | 1.58 | 3.87 | 5.55  |
| ACA1    | 1.14 | 1.38 | 1.18 | 1.08 | 1.09 | 2.03 | 1.87 | 4.08 | 4.23  |
| YOR084W | 1.01 | 1.49 | 1.25 | 1.00 | 1.35 | 2.11 | 2.24 | 3.40 | 15.54 |
| SHO1    | 1.39 | 2.60 | 1.72 | 0.81 | 0.67 | 1.75 | 0.96 | 3.91 | 8.85  |
| GNP1    | 1.13 | 1.55 | 1.67 | 0.94 | 0.99 | 1.76 | 2.72 | 3.05 | 4.49  |
| YNL081C | 0.99 | 0.92 | 0.88 | 0.55 | 0.58 | 3.08 | 2.35 | 4.45 | 12.04 |
| YKL161C | 1.01 | 1.56 | 1.56 | 1.23 | 1.18 | 2.24 | 2.16 | 2.85 | 6.02  |
| ISY1    | 1.05 | 0.97 | 1.21 | 1.13 | 1.33 | 2.01 | 2.89 | 3.17 | 9.46  |
| BEM2    | 0.67 | 1.14 | 0.49 | 0.73 | 1.12 | 2.07 | 2.34 | 5.19 | 14.32 |
| PCP1    | 1.05 | 1.27 | 1.24 | 1.40 | 1.71 | 1.51 | 2.04 | 3.51 | 5.83  |
| CTF19   | 1.25 | 1.25 | 0.99 | 1.01 | 1.01 | 1.36 | 3.54 | 3.28 | 6.82  |
| ACF2    | 1.17 | 1.33 | 1.49 | 1.63 | 1.43 | 1.59 | 1.38 | 3.66 | 4.66  |
| MCM21   | 0.81 | 1.40 | 1.28 | 0.88 | 1.28 | 1.67 | 2.12 | 4.16 | 13.63 |
| BUD13   | 0.97 | 1.34 | 0.45 | 0.93 | 1.25 | 3.07 | 2.35 | 3.24 | 11.70 |
| YPR097W | 0.98 | 0.85 | 0.98 | 0.87 | 1.01 | 1.95 | 4.35 | 2.60 | 4.74  |
| YCR025C | 1.05 | 3.71 | 0.81 | 1.05 | 2.82 | 0.81 | 1.05 | 2.27 | 3.81  |
| MDM31   | 0.74 | 1.14 | 0.64 | 1.02 | 1.31 | 1.35 | 1.94 | 5.41 | 6.32  |

|           |      |      |      |      |      |      |      |      |       |
|-----------|------|------|------|------|------|------|------|------|-------|
| TOM5      | 1.68 | 0.60 | 0.62 | 0.93 | 1.65 | 1.13 | 3.04 | 3.84 | 12.12 |
| IBD2      | 0.76 | 1.16 | 0.72 | 0.67 | 0.80 | 2.51 | 2.34 | 4.53 | 26.35 |
| YHR132W-A | 0.95 | 0.71 | 1.09 | 1.18 | 0.97 | 1.45 | 1.57 | 5.55 | 11.24 |
| PHO2      | 1.02 | 1.15 | 1.06 | 1.01 | 1.11 | 1.86 | 3.90 | 2.33 | 4.45  |
| SHE4      | 0.81 | 1.25 | 0.89 | 0.85 | 1.58 | 2.35 | 1.89 | 3.61 | 9.21  |
| JIP3      | 0.77 | 1.63 | 0.89 | 1.20 | 1.16 | 1.97 | 3.05 | 2.51 | 3.72  |
| CLB4      | 0.98 | 0.98 | 1.01 | 1.03 | 1.13 | 1.17 | 2.20 | 4.63 | 19.30 |
| YDR010C   | 0.84 | 0.96 | 0.67 | 0.81 | 0.69 | 1.57 | 3.10 | 4.44 | 8.40  |
| YPR045C   | 1.17 | 1.01 | 0.87 | 0.92 | 0.83 | 1.25 | 1.50 | 5.50 | 10.78 |
| HST1      | 1.30 | 0.96 | 1.09 | 1.15 | 1.57 | 1.25 | 1.91 | 3.73 | 3.72  |
| DPH2      | 1.20 | 0.94 | 1.27 | 0.81 | 0.96 | 1.48 | 1.76 | 4.50 | 5.77  |
| PRY2      | 0.85 | 1.52 | 1.15 | 1.32 | 1.66 | 1.77 | 1.53 | 2.97 | 4.95  |
| LIP5      | 0.75 | 1.50 | 1.10 | 0.30 | 1.42 | 1.96 | 1.10 | 4.51 | 13.22 |
| SSP1      | 1.20 | 1.30 | 1.16 | 1.26 | 1.39 | 1.65 | 2.09 | 2.58 | 5.82  |
| YDR133C   | 0.57 | 1.29 | 1.36 | 1.23 | 1.23 | 2.44 | 1.28 | 3.20 | 7.46  |
| MCM16     | 0.97 | 1.17 | 0.90 | 0.86 | 1.24 | 1.45 | 1.71 | 4.22 | 9.91  |
| YBR042C   | 1.28 | 1.46 | 1.28 | 1.24 | 1.49 | 1.35 | 1.48 | 2.91 | 4.24  |
| YMR041C   | 0.94 | 1.16 | 1.23 | 0.93 | 0.91 | 2.24 | 2.34 | 2.75 | 4.74  |
| MAK3      | 0.99 | 1.17 | 1.46 | 1.09 | 1.40 | 1.35 | 2.34 | 2.67 | 3.66  |
| MEU1      | 0.91 | 1.77 | 0.93 | 0.91 | 1.20 | 1.65 | 1.61 | 3.41 | 5.58  |
| ALK1      | 0.93 | 1.15 | 0.93 | 1.11 | 1.09 | 1.76 | 1.60 | 3.81 | 4.57  |
| YJL046W   | 0.76 | 0.79 | 1.13 | 0.61 | 0.92 | 2.26 | 2.74 | 2.93 | 2.63  |
| MBF1      | 0.73 | 1.28 | 1.11 | 0.79 | 1.12 | 1.64 | 2.95 | 2.49 | 2.44  |
| TES1      | 1.12 | 1.35 | 1.42 | 0.93 | 1.09 | 2.30 | 1.34 | 2.49 | 4.82  |
| MRS2      | 0.74 | 1.08 | 0.60 | 0.71 | 0.77 | 1.92 | 2.61 | 3.60 | 5.41  |
| YEL014C   | 0.97 | 0.85 | 1.25 | 1.07 | 0.87 | 1.45 | 1.47 | 4.08 | 3.93  |
| PBP2      | 1.10 | 0.64 | 0.77 | 0.77 | 0.86 | 0.80 | 1.06 | 5.99 | 4.77  |
| YPR076W   | 0.77 | 1.26 | 1.64 | 0.70 | 0.96 | 1.64 | 1.61 | 3.41 | 3.65  |
| YDR219C   | 1.08 | 0.63 | 0.99 | 1.05 | 1.38 | 1.78 | 2.08 | 2.74 | 6.57  |
| BMH1      | 0.92 | 1.28 | 1.30 | 0.88 | 0.95 | 1.18 | 1.53 | 3.66 | 11.63 |
| CEM1      | 0.83 | 0.72 | 0.54 | 0.78 | 0.73 | 1.86 | 2.17 | 3.99 | 7.81  |
| YBR269C   | 0.64 | 0.93 | 1.02 | 1.15 | 0.93 | 2.10 | 2.45 | 2.35 | 4.01  |
| RRD2      | 0.82 | 0.73 | 0.55 | 0.66 | 0.90 | 0.99 | 2.24 | 4.63 | 4.07  |
| YDR134C   | 0.83 | 1.20 | 0.96 | 1.16 | 1.37 | 1.52 | 2.07 | 2.34 | 7.58  |
| OM45      | 1.09 | 1.07 | 1.25 | 0.92 | 1.22 | 1.27 | 2.05 | 2.54 | 5.04  |
| YMR160W   | 0.95 | 1.34 | 1.07 | 0.96 | 1.28 | 1.40 | 1.62 | 2.74 | 4.26  |
| PPT2      | 1.14 | 1.06 | 1.38 | 0.83 | 0.94 | 1.49 | 1.90 | 2.60 | 3.64  |
| YPL017C   | 1.36 | 1.12 | 0.87 | 0.83 | 1.13 | 1.03 | 1.87 | 3.07 | 6.40  |
| INP54     | 0.54 | 1.32 | 1.62 | 1.00 | 0.93 | 1.36 | 1.92 | 2.61 | 6.83  |
| FIL1      | 1.08 | 0.81 | 0.60 | 0.67 | 1.03 | 1.29 | 1.37 | 4.35 | 11.01 |
| ASP1      | 0.99 | 1.05 | 1.48 | 1.03 | 0.96 | 1.80 | 1.60 | 2.26 | 5.31  |
| EAF5      | 1.10 | 1.09 | 0.95 | 0.84 | 0.66 | 1.54 | 1.24 | 3.71 | 10.68 |
| PNP1      | 1.07 | 0.79 | 0.84 | 0.92 | 0.91 | 1.59 | 2.08 | 2.92 | 6.37  |
| LCB5      | 1.34 | 1.05 | 1.23 | 1.06 | 1.23 | 1.38 | 1.44 | 2.38 | 6.27  |

|         |      |      |      |      |      |      |      |      |       |
|---------|------|------|------|------|------|------|------|------|-------|
| YIL028W | 1.35 | 1.28 | 0.92 | 0.88 | 1.09 | 1.61 | 1.37 | 2.61 | 4.04  |
| BRR1    | 0.95 | 1.03 | 0.52 | 0.64 | 0.63 | 1.24 | 1.95 | 3.98 | 6.87  |
| TRM1    | 1.06 | 1.88 | 0.97 | 0.84 | 0.75 | 1.21 | 1.41 | 2.80 | 8.33  |
| ISW2    | 1.12 | 0.98 | 0.72 | 0.95 | 0.62 | 2.05 | 1.02 | 3.34 | 3.99  |
| RPN10   | 0.82 | 0.66 | 0.63 | 0.79 | 0.93 | 0.52 | 0.56 | 5.88 | 7.40  |
| KSS1    | 1.01 | 1.23 | 1.15 | 0.90 | 0.96 | 0.86 | 1.08 | 3.60 | 5.78  |
| YKR032W | 0.85 | 1.90 | 0.99 | 1.11 | 1.20 | 1.25 | 1.22 | 2.25 | 4.56  |
| CTF8    | 0.85 | 1.24 | 1.07 | 1.01 | 1.03 | 1.12 | 1.36 | 3.08 | 5.56  |
| CSN12   | 0.86 | 1.86 | 0.80 | 1.02 | 0.82 | 1.02 | 1.57 | 2.75 | 4.37  |
| SYF2    | 0.81 | 1.09 | 1.33 | 1.04 | 1.42 | 1.55 | 1.17 | 2.22 | 3.65  |
| DOG2    | 0.89 | 1.22 | 1.20 | 0.94 | 0.94 | 1.17 | 1.57 | 2.58 | 3.22  |
| ELM1    | 1.00 | 0.84 | 0.71 | 0.76 | 0.72 | 1.22 | 1.33 | 3.35 | 18.98 |
| LEA1    | 1.10 | 0.57 | 0.36 | 0.50 | 0.43 | 1.26 | 1.25 | 4.28 | 4.87  |
| DSE4    | 0.91 | 0.90 | 0.90 | 1.16 | 0.97 | 0.86 | 1.58 | 2.35 | 4.04  |
| YCK3    | 1.07 | 0.45 | 0.48 | 0.62 | 0.68 | 1.14 | 2.03 | 2.88 | 4.50  |
| SPO1    | 0.77 | 1.27 | 1.13 | 0.90 | 0.60 | 0.92 | 0.77 | 2.94 | 4.22  |

Notes: Benomyl-sensitivity of each freshly generated *MATa* haploid YKO mutant was measured by microarray analysis of the representing TAGs in the experimental pool (treated with drug of the indicated concentrations) as compared to a control (drug-untreated) pool. A high C/E ratio indicates that the mutant is benomyl-sensitive and was shaded in gray.

**Table S3. Heterozygous diploid YKO strains excluded from the pool**

| <b>ORF name</b> | <b>plate</b> | <b>His<sup>+</sup><sup>a</sup></b> | <b>Mates as<br/><i>MATa</i><sup>b</sup></b> | <b>Papillates as<br/><i>MATa</i><sup>b</sup></b> | <b>Mates as<br/><i>MATa</i><sup>b</sup></b> | <b>Papillates as<br/><i>MATa</i><sup>b</sup></b> |
|-----------------|--------------|------------------------------------|---|--|---|--|
| YAL051W         | 201          | 1                                  | 0   | 0  | 0   | 0  |
| YAL035W         | 201          | 1                                  | 0   | 0  | 0   | 0  |
| YAL017W         | 201          | 1                                  | 0   | 0  | 0   | 0  |
| YAL005C         | 201          | 0                                  | 1   | 0  | 0   | 0  |
| YAR002W         | 201          | 0                                  | 0   | 0  | 1   | 0  |
| YAR018C         | 201          | 0                                  | 1   | 0  | 0   | 0  |
| YML089C         | 203          | 1                                  | 0   | 0  | 0   | 0  |
| YML088W         | 203          | 1                                  | 0   | 0  | 0   | 0  |
| YML087C         | 203          | 1                                  | 0   | 0  | 0   | 0  |
| YML086C         | 203          | 1                                  | 0   | 0  | 0   | 0  |
| YML084W         | 203          | 1                                  | 0   | 0  | 0   | 0  |
| YML083C         | 203          | 1                                  | 0   | 0  | 0   | 0  |
| YML082W         | 203          | 1                                  | 0   | 0  | 0   | 0  |
| YML081W         | 203          | 1                                  | 0   | 0  | 0   | 0  |
| YML080W         | 203          | 1                                  | 0   | 0  | 0   | 0  |
| YML079W         | 203          | 1                                  | 0   | 0  | 0   | 0  |
| YML078W         | 203          | 1                                  | 0   | 0  | 0   | 0  |
| YMR223W         | 205          | 1                                  | 0   | 0  | 0   | 0  |
| YPL273W         | 209          | 0                                  | 1   | 0  | 0   | 0  |
| YPL202C         | 210          | 0                                  | 0   | 0  | 1   | 0  |
| YER068W         | 215          | 0                                  | 0   | 1  | 0   | 0  |
| YHR080C         | 217          | 0                                  | 0   | 0  | 0   | 1  |
| YHR094C         | 217          | 0                                  | 0   | 0  | 0   | 1  |
| YCL014W         | 218          | 0                                  | 0   | 0  | 0   | 1  |
| YCL044C         | 219          | 0                                  | 0   | 1  | 0   | 0  |
| YLR151C         | 219          | 0                                  | 0   | 0  | 0   | 1  |
| YGR051C         | 221          | 0                                  | 0   | 0  | 0   | 1  |
| YGR065C         | 221          | 0                                  | 0   | 0  | 0   | 1  |
| YOR138C         | 223          | 0                                  | 0   | 1  | 0   | 0  |
| YOR143C         | 223          | 0                                  | 0   | 0  | 0   | 1  |
| YGL215W         | 228          | 0                                  | 0   | 0  | 1   | 0  |
| YIL046W         | 237          | 0                                  | 0   | 1  | 0   | 0  |
| YFL030W         | 237          | 0                                  | 0   | 0  | 0   | 1  |
| YJL111W         | 242          | 0                                  | 1   | 0  | 0   | 1  |
| YJR148W         | 246          | 0                                  | 0   | 0  | 0   | 1  |
| YLR343W         | 246          | 0                                  | 0   | 0  | 0   | 1  |
| YML031W         | 246          | 0                                  | 0   | 0  | 0   | 1  |
| YML046W         | 246          | 0                                  | 0   | 0  | 0   | 1  |
| YML067C         | 246          | 0                                  | 0   | 0  | 0   | 1  |
| YDL217C         | 246          | 0                                  | 0   | 1  | 0   | 0  |



|         |     |   |   |   |   |   |
|---------|-----|---|---|---|---|---|
| YDR004W | 246 | 0 | 0 | 1 | 0 | 0 |
| YDR006C | 246 | 0 | 0 | 0 | 0 | 1 |
| YNR017W | 252 | 0 | 0 | 0 | 0 | 1 |
| YBR159W | 258 | 0 | 0 | 1 | 0 | 0 |
| YBR161W | 258 | 0 | 0 | 0 | 0 | 1 |
| YIL102C | 262 | 0 | 0 | 0 | 1 | 0 |
| YIL051C | 262 | 0 | 1 | 0 | 0 | 0 |
| YJL166W | 280 | 0 | 1 | 0 | 0 | 0 |
| YLR069C | 280 | 0 | 0 | 0 | 0 | 1 |
| YML060W | 280 | 0 | 0 | 0 | 0 | 1 |
| YML052W | 280 | 0 | 0 | 0 | 0 | 1 |
| YMR266W | 280 | 0 | 1 | 0 | 0 | 0 |
| YNL331C | 280 | 0 | 1 | 0 | 0 | 0 |
| YNL262W | 280 | 0 | 0 | 0 | 1 | 0 |
| YNL251C | 280 | 0 | 0 | 0 | 0 | 0 |
| YOR014W | 280 | 0 | 0 | 0 | 0 | 1 |
| YOR078W | 280 | 0 | 0 | 0 | 0 | 1 |
| YOR089C | 280 | 0 | 0 | 1 | 0 | 0 |
| YOR093C | 280 | 0 | 0 | 1 | 0 | 0 |
| YOR095C | 280 | 0 | 0 | 1 | 0 | 1 |
| YOR097C | 280 | 0 | 0 | 1 | 0 | 1 |
| YOR106W | 280 | 0 | 0 | 1 | 0 | 1 |
| YOR132W | 280 | 0 | 0 | 1 | 0 | 0 |
| YOR149C | 280 | 0 | 0 | 1 | 0 | 0 |
| YOR202W | 280 | 0 | 0 | 1 | 0 | 0 |
| YOR204W | 280 | 0 | 0 | 1 | 0 | 1 |
| YOR278W | 280 | 0 | 1 | 0 | 0 | 1 |
| YOR279C | 280 | 0 | 0 | 1 | 0 | 1 |
| YOR281C | 280 | 0 | 0 | 1 | 0 | 1 |
| YDR077W | 280 | 0 | 1 | 0 | 0 | 0 |
| YEL065W | 280 | 0 | 0 | 1 | 0 | 1 |
| YER007W | 280 | 0 | 0 | 0 | 0 | 0 |
| YER058W | 280 | 0 | 0 | 0 | 0 | 0 |
| YER086W | 280 | 0 | 0 | 0 | 0 | 0 |
| YGR096W | 280 | 0 | 0 | 0 | 0 | 0 |
| YGR111W | 280 | 0 | 0 | 0 | 0 | 0 |
| YGR186W | 280 | 0 | 0 | 0 | 0 | 0 |
| YGR198W | 280 | 0 | 0 | 0 | 0 | 0 |
| YGR199W | 280 | 0 | 0 | 0 | 0 | 0 |
| YGR206W | 280 | 0 | 0 | 0 | 0 | 0 |
| YHR148W | 280 | 0 | 0 | 0 | 0 | 0 |
| YHR160C | 280 | 0 | 0 | 0 | 0 | 0 |
| YHR195W | 280 | 0 | 0 | 0 | 0 | 0 |
| YHR207C | 280 | 0 | 0 | 1 | 0 | 1 |

|         |     |   |   |   |   |   |
|---------|-----|---|---|---|---|---|
| YBL023C | 281 | 0 | 1 | 0 | 0 | 1 |
| YBL091C | 281 | 0 | 1 | 0 | 0 | 0 |
| YBR023C | 281 | 0 | 0 | 1 | 0 | 1 |
| YFL008W | 281 | 0 | 1 | 0 | 0 | 0 |
| YER157W | 281 | 0 | 1 | 0 | 0 | 0 |
| YMR075W | 281 | 0 | 1 | 0 | 0 | 1 |
| YNL004W | 281 | 0 | 0 | 0 | 1 | 0 |
| YNL169C | 281 | 0 | 0 | 0 | 1 | 0 |
| YIL102C | 281 | 0 | 0 | 1 | 0 | 1 |
| YOR162C | 281 | 0 | 0 | 1 | 0 | 1 |
| YBL021C | 281 | 0 | 0 | 0 | 1 | 0 |
| YDR501W | 281 | 0 | 0 | 0 | 1 | 0 |
| YKR078W | 281 | 0 | 0 | 1 | 1 | 0 |
| YBR272C | 281 | 0 | 0 | 1 | 0 | 1 |
| YFL039C | 281 | 0 | 1 | 0 | 1 | 0 |
| YCL026C | 281 | 0 | 0 | 1 | 0 | 1 |

Notes: “1” means “yes”, whereas “0” means “no” as results of strain quality control. “a” His<sup>+</sup> was defined by the growth of the strain in synthetic medium lacking L-Histidine. “b” mating phenotypes were tested with wild-type strains BY4741a, BY4742□, and BY4743a/□ as control. His<sup>+</sup> strains and mutants behaved like haploid in mating type tests were excluded from the pool used in this study. All strains on plates 280 and 281 were excluded because most strains on these two plates failed in the quality tests. In total, 100 of the 5,996 strains were not included in the heterozygous diploid mutant pool used in this study. The identity of all other 5,896 strains used in this study can be obtained from Research genetics.

**Table S4. Internal controls of the *cin8* dSLAM screens**

| ORF     | Gene Name | Gene description         | single/double hybridization |          |          |            |
|---------|-----------|--------------------------|-----------------------------|----------|----------|------------|
|         |           |                          | ratio #1                    | ratio #2 | ratio #3 | mean ratio |
| YKL216W | URA1      | Uracil synthesis         | 22.17                       | 39.85    | 22.69    | 28.24      |
| YJL130C | URA2      | Uracil synthesis         | 10.68                       | 6.27     | 6.35     | 7.77       |
| YLR420W | URA4      | Uracil synthesis         | 132.08                      | 45.14    | 141.79   | 106.34     |
| YML106W | URA5      | Uracil synthesis         | 27.89                       | 44.94    | 51.26    | 41.36      |
| YLR014C | PPR1      | Uracil synthesis         | 6.98                        | 7.87     | 4.95     | 6.60       |
| YEL060C | PRB1      | right to <i>CIN8</i> ORF | 10.52                       | 8.68     | 9.10     | 9.43       |
| YEL061C | CIN8      | query gene               | 0.19                        | 0.20     | 0.16     | 0.18       |
| YEL062W | NPR2      | left to <i>CIN8</i> ORF  | 13.73                       | 12.81    | 7.81     | 11.45      |

Notes: In each *cin8* dSLAM experiment, useful UPTAG and DOWNTAG ratios were averaged to obtain a single/double (C/E) hybridization ratio for each YKO. A mean ratio was taken among all three experiments.

**Table S5. A summary of the *cin8* dSLAM screen.**

| ORF       | Gene Name | single/double hybridization |          |          |            | confirmation |            |
|-----------|-----------|-----------------------------|----------|----------|------------|--------------|------------|
|           |           | ratio #1                    | ratio #2 | ratio #3 | mean ratio | RSA          | Tetrad     |
| YJL013C   | MAD3      | 52.95                       | 120.19   | 37.11    | 70.08      | SL           | SL         |
| YER177W   | BMH1      | 5.49                        | 95.83    | 9.50     | 36.94      | SL           | SL         |
| YBL063W   | KIP1      | 46.98                       | 31.25    | 20.40    | 32.88      | SL           | SL         |
| YPR119W   | CLB2      | 20.58                       | 20.23    | 31.28    | 24.03      | SL           | SL         |
| YCR076C   | □         | 3.27                        | 61.16    | 5.27     | 23.23      | NO           | undecisive |
| YER135C   | □         | 31.68                       | 20.60    | 17.22    | 23.17      | SL           | undecisive |
| YMR055C   | BUB2      | 23.63                       | 15.68    | 20.15    | 19.82      | SL           | SL         |
| YHR013C   | ARD1      | 22.28                       | 18.73    | 14.12    | 18.38      | SL           | SL         |
| YEL003W   | GIM4      | 21.06                       | 19.10    | 13.81    | 17.99      | SF           | SL         |
| YOR073W   | SGO1      | 10.62                       | 10.47    | 31.89    | 17.66      | SL           | SL         |
| YDR149C   | □         | 21.66                       | 21.16    | 9.67     | 17.50      | SL           | SL         |
| YJL179W   | PFD1      | 13.61                       | 12.27    | 14.54    | 13.47      | SL           | SL         |
| YMR299C   | DYN3      | 8.12                        | 16.52    | 11.67    | 12.10      | SL           | SL         |
| YMR294W   | JNM1      | 15.54                       | 11.12    | 7.41     | 11.36      | SL           | SL         |
| YMR078C   | CTF18     | 11.35                       | 12.24    | 10.35    | 11.32      | SL           | SL         |
| YLR226W   | BUR2      | 16.13                       | 5.39     | 12.41    | 11.31      | SL           | SL         |
| YNL016W   | PUB1      | 3.99                        | 2.97     | 25.93    | 10.96      | SF           | SF         |
| YKR061W   | KTR2      | 8.67                        | 5.80     | 15.83    | 10.10      | SF           | SL         |
| YDR488C   | PAC11     | 5.66                        | 13.34    | 11.00    | 10.00      | SL           | SL         |
| YPL174C   | NIP100    | 10.44                       | 7.93     | 10.55    | 9.64       | SL           | SL         |
| YHR191C   | CTF8      | 8.61                        | 11.27    | 8.58     | 9.49       | SL           | SL         |
| YML010W-A | □         | 3.06                        | 7.86     | 17.45    | 9.46       | SL           | SL         |
| YDR150W   | NUM1      | 7.85                        | 10.51    | 9.05     | 9.14       | SL           | SL         |
| YOR269W   | PAC1      | 10.41                       | 7.53     | 7.90     | 8.61       | SL           | SL         |
| YBR194W   | SOY1      | 3.52                        | 10.73    | 11.56    | 8.60       | SF           | SF         |
| YLR418C   | CDC73     | 8.16                        | 9.40     | 7.50     | 8.35       | SL           | SL         |
| YDR432W   | NPL3      | 16.21                       | 3.81     | 3.15     | 7.72       | SL           | SL         |
| YLL049W   | □         | 8.96                        | 7.84     | 5.38     | 7.39       | SL           | SL         |
| YCR086W   | CSM1      | 7.47                        | 10.16    | 4.23     | 7.29       | SL           | SL         |
| YDR435C   | PPM1      | 7.55                        | 7.05     | 6.68     | 7.09       | SF           | SF         |
| YLR210W   | CLB4      | 6.69                        | 7.92     | 6.41     | 7.01       | SL           | SL         |
| YBL058W   | SHP1      | 6.30                        | 10.21    | 4.27     | 6.93       | SL           | SL         |
| YNL147W   | LSM7      | 6.54                        | 5.87     | 8.35     | 6.92       | SL           | SL         |
| YOR195W   | SLK19     | 7.01                        | 6.44     | 7.11     | 6.85       | SL           | SL         |
| YPL008W   | CHL1      | 7.16                        | 6.63     | 5.81     | 6.53       | SL           | SL         |
| YBR140C   | IRA1      | 7.86                        | 6.42     | 5.19     | 6.49       | SF           | SL         |
| YCL016C   | DCC1      | 5.22                        | 8.20     | 5.86     | 6.43       | SL           | SL         |
| YGL086W   | MAD1      | 7.16                        | 5.52     | 6.43     | 6.37       | SL           | SL         |
| YGL003C   | CDH1      | 4.64                        | 5.09     | 9.20     | 6.31       | SL           | SL         |

|         |        |      |      |      |      |    |            |
|---------|--------|------|------|------|------|----|------------|
| YGR188C | BUB1   | 5.91 | 4.12 | 8.88 | 6.30 | SL | SL         |
| YJL030W | MAD2   | 6.95 | 5.02 | 6.78 | 6.25 | SL | SL         |
| YOR349W | CIN1   | 5.00 | 6.99 | 6.21 | 6.07 | SL | SL         |
| YJL124C | LSM1   | 5.79 | 5.43 | 6.55 | 5.92 | SL | SL         |
| YNL153C | GIM3   | 6.87 | 4.50 | 5.59 | 5.65 | SL | SL         |
| YML094W | GIM5   | 5.71 | 6.51 | 4.60 | 5.61 | SL | SL         |
| YKL057C | NUP120 | 6.82 | 4.30 | 5.57 | 5.56 | SL | SL         |
| YNR052C | POP2   | 6.15 | 4.98 | 5.10 | 5.41 | SL | SL         |
| YIL125W | KGD1   | 6.59 | 4.40 | 5.17 | 5.39 | NO | NO         |
| YGL031C | RPL24A | 9.72 | 2.68 | 3.55 | 5.32 | NO | NO         |
| YAL024C | LTE1   | 5.01 | 4.93 | 5.94 | 5.29 | SL | SL         |
| YHR041C | SRB2   | 3.48 | 5.48 | 6.72 | 5.22 | SL | SL         |
| YJR053W | BFA1   | 4.04 | 5.30 | 6.33 | 5.22 | SL | SL         |
| YBR171W | SEC66  | 7.43 | 4.14 | 3.88 | 5.15 | SF | SL         |
| YDL040C | NAT1   | 4.97 | 5.58 | 4.76 | 5.10 | SL | SL         |
| YKR010C | TOF2   | 4.85 | 5.58 | 4.83 | 5.09 | SF | SF         |
| YHL019C | APM2   | 5.20 | 5.13 | 4.76 | 5.03 | NO | NO         |
| YMR116C | ASC1   | 6.63 | 5.00 | 2.83 | 4.82 | NO | NO         |
| YDL134C | PPH21  | 5.47 | 4.90 | 4.03 | 4.80 | SF | SF         |
| YDR424C | DYN2   | 5.91 | 3.83 | 4.37 | 4.70 | SF | SL         |
| YGR078C | PAC10  | 6.17 | 4.45 | 3.45 | 4.69 | SL | SL         |
| YHR129C | ARP1   | 4.83 | 3.73 | 5.49 | 4.68 | SF | SL         |
| YLR268W | SEC22  | 4.97 | 4.85 | 4.11 | 4.64 | SL | SL         |
| YNL269W | □      | 4.38 | 4.73 | 4.75 | 4.62 | NO | undecisive |
| YDL229W | SFB1   | 3.39 | 4.85 | 5.45 | 4.56 | NO | NO         |
| YKL191W | DPH2   | 2.65 | 7.85 | 2.99 | 4.50 | NO | NO         |
| YLR200W | YKE2   | 4.79 | 4.70 | 3.84 | 4.44 | SF | SL         |
| YNL171C | □      | 7.14 | 3.71 | 2.31 | 4.39 | SL | SL         |
| YNL225C | CNM67  | 7.32 | 2.66 | 3.07 | 4.35 | SL | SL         |
| YLR338W | □      | 4.71 | 4.56 | 3.56 | 4.28 | SF | SF         |
| YER019W | ISC1   | 3.36 | 4.57 | 4.78 | 4.24 | SF | SF         |
| YIL137C | □      | 3.59 | 4.34 | 4.67 | 4.20 | NO | NO         |
| YBR282W | MRPL27 | 6.32 | 4.02 | 2.25 | 4.20 | NO | NO         |
| YOR360C | PDE2   | 4.85 | 4.05 | 3.65 | 4.18 | SF | SF         |
| YKL048C | ELM1   | 3.95 | 4.03 | 4.51 | 4.17 | SF | SL         |
| YOR058C | ASE1   | 3.63 | 4.02 | 4.83 | 4.16 | SL | SL         |
| YEL044W | IES6   | 5.32 | 3.86 | 3.22 | 4.14 | SF | SF         |
| YNL273W | TOF1   | 4.33 | 4.19 | 3.41 | 3.98 | SF | SL         |
| YKL176C | LST4   | 5.91 | 3.18 | 2.80 | 3.96 | SF | SF         |
| YGL167C | PMR1   | 3.53 | 3.89 | 4.26 | 3.89 | SF | SL         |
| YLR315W | NKP2   | 4.15 | 3.26 | 4.26 | 3.89 | SF | SF         |
| YLR254C | □      | 3.66 | 4.14 | 3.82 | 3.87 | SF | SL         |
| YIL114C | POR2   | 4.52 | 3.33 | 3.65 | 3.83 | NO | NO         |

|           |        |      |      |      |      |    |    |
|-----------|--------|------|------|------|------|----|----|
| YIL039W   | □      | 2.56 | 3.43 | 5.32 | 3.77 | NO | NO |
| YDL234C   | GYP7   | 4.61 | 2.92 | 3.78 | 3.77 | NO | NO |
| YPR135W   | CTF4   | 3.56 | 3.34 | 4.35 | 3.75 | SL | SL |
| YGR165W   | MRPS35 | 3.27 | 4.96 | 3.00 | 3.74 | NO | NO |
| YIL153W   | RRD1   | 5.53 | 3.21 | 2.39 | 3.71 | SF | SL |
| YGR285C   | ZUO1   | 3.33 | 3.34 | 4.39 | 3.69 | SF | SL |
| YOR026W   | BUB3   | 4.18 | 2.58 | 4.28 | 3.68 | SL | SL |
| YFR055W   | □      | 2.27 | 3.03 | 5.70 | 3.67 | NO | □  |
| YGL219C   | MDM34  | 5.01 | 3.53 | 2.17 | 3.57 | SF | SF |
| YGL153W   | PEX14  | 2.58 | 4.14 | 3.89 | 3.54 | SF | NO |
| YPL241C   | CIN2   | 3.03 | 4.07 | 3.33 | 3.48 | SL | SL |
| YNL307C   | MCK1   | 2.69 | 3.54 | 4.00 | 3.41 | SF | SL |
| YDR532C   | □      | 3.26 | 2.53 | 4.39 | 3.40 | SL | □  |
| YGR104C   | SRB5   | 4.31 | 3.40 | 2.44 | 3.38 | SF | SL |
| YIL165C   | □      | 3.31 | 4.10 | 2.72 | 3.38 | NO | □  |
| YEL017W   | GTT3   | 2.45 | 5.41 | 2.25 | 3.37 | NO | □  |
| YKL204W   | EAP1   | 3.60 | 2.42 | 3.93 | 3.32 | SF | SF |
| YML016C   | PPZ1   | 3.42 | 2.84 | 3.62 | 3.29 | SF | SF |
| YKR047W   | □      | 3.63 | 3.30 | 2.94 | 3.29 | NO | □  |
| YKR059W   | TIF1   | 2.79 | 3.70 | 3.31 | 3.27 | SF | SF |
| YMR048W   | CSM3   | 2.64 | 3.30 | 3.79 | 3.24 | SL | SL |
| YOL076W   | MDM20  | 3.45 | 3.25 | 3.00 | 3.23 | SF | SF |
| YIL055C   | □      | 3.02 | 3.53 | 3.12 | 3.23 | NO | □  |
| YPR141C   | KAR3   | 2.48 | 4.33 | 2.85 | 3.22 | SF | SL |
| YIL154C   | IMP2   | 3.49 | 3.37 | 2.76 | 3.21 | NO | □  |
| YJL169W   | □      | 3.16 | 3.35 | 3.08 | 3.20 | SF | SF |
| YOL081W   | IRA2   | 4.18 | 2.89 | 2.46 | 3.17 | SF | SF |
| YGR238C   | KEL2   | 3.03 | 2.87 | 3.62 | 3.17 | NO | □  |
| YIL141W   | □      | 2.88 | 3.93 | 2.64 | 3.15 | NO | □  |
| YPL018W   | CTF19  | 3.16 | 3.01 | 3.23 | 3.14 | SL | SL |
| YJR102C   | VPS25  | 3.62 | 2.85 | 2.88 | 3.12 | SF | SF |
| YHR064C   | SFZ1   | 2.79 | 2.98 | 3.53 | 3.10 | SF | SF |
| YLR417W   | VPS36  | 2.91 | 2.97 | 3.32 | 3.07 | NO | □  |
| YGR229C   | SMI1   | 3.80 | 3.02 | 2.14 | 2.98 | SF | SL |
| YIL073C   | SPO22  | 2.81 | 3.07 | 3.05 | 2.97 | NO | □  |
| YBR258C   | SHG1   | 2.91 | 3.54 | 2.46 | 2.97 | NO | □  |
| YKR054C   | DYN1   | 2.86 | 3.23 | 2.80 | 2.97 | SF | SL |
| YIL160C   | POT1   | 2.59 | 3.00 | 3.28 | 2.95 | NO | □  |
| YML094C-A | □      | 3.62 | 2.83 | 2.34 | 2.93 | NO | SL |
| YNL332W   | THI12  | 2.66 | 2.77 | 3.34 | 2.92 | NO | □  |
| YLR388W   | RPS29A | 3.57 | 2.45 | 2.72 | 2.91 | SF | SF |
| YOR288C   | MPD1   | 2.82 | 3.11 | 2.78 | 2.90 | NO | □  |
| YIL159W   | BNR1   | 2.93 | 2.63 | 3.13 | 2.90 | NO | □  |

|         |        |      |      |      |      |    |    |
|---------|--------|------|------|------|------|----|----|
| YGL127C | SOH1   | 2.22 | 3.91 | 2.45 | 2.86 | NO | □  |
| YLL033W | □      | 3.14 | 2.58 | 2.85 | 2.85 | NO | □  |
| YGL216W | KIP3   | 2.41 | 3.08 | 3.06 | 2.85 | SF | SL |
| YER095W | RAD51  | 2.01 | 3.27 | 3.27 | 2.85 | NO | NO |
| YNL136W | □      | 2.84 | 2.51 | 3.17 | 2.84 | SF | SF |
| YDR334W | SWR1   | 2.51 | 3.28 | 2.61 | 2.80 | SF | SF |
| YNL298W | CLA4   | 2.92 | 3.20 | 2.25 | 2.79 | SF | SL |
| YLR084C | RAX2   | 3.09 | 2.20 | 3.05 | 2.78 | NO | □  |
| YGR237C | □      | 2.71 | 2.85 | 2.67 | 2.74 | NO | □  |
| YIL097W | FYV10  | 2.40 | 2.78 | 2.93 | 2.70 | NO | □  |
| YGR077C | PEX8   | 2.65 | 2.64 | 2.83 | 2.70 | SF | NO |
| YOL012C | HTZ1   | 2.44 | 2.43 | 3.20 | 2.69 | SF | SF |
| YJL168C | SET2   | 2.82 | 2.92 | 2.33 | 2.69 | SF | SF |
| YFR039C | □      | 2.58 | 2.82 | 2.56 | 2.66 | NO | □  |
| YHR200W | RPN10  | 2.22 | 2.87 | 2.86 | 2.65 | NO | □  |
| YER040W | GLN3   | 2.41 | 2.76 | 2.77 | 2.65 | NO | □  |
| YOR358W | HAP5   | 2.24 | 2.85 | 2.84 | 2.64 | NO | □  |
| YOR243C | PUS7   | 2.78 | 2.67 | 2.46 | 2.64 | SF | SF |
| YOL159C | □      | 2.52 | 2.77 | 2.62 | 2.64 | NO | □  |
| YDL115C | IWR1   | 2.24 | 3.52 | 2.14 | 2.63 | SF | SL |
| YHL023C | RMD11  | 2.72 | 3.11 | 2.05 | 2.63 | NO | □  |
| YPL079W | RPL21B | 2.75 | 2.98 | 2.15 | 2.62 | NO | □  |
| YNL314W | DAL82  | 2.60 | 2.54 | 2.72 | 2.62 | NO | □  |
| YAL010C | MDM10  | 2.86 | 2.53 | 2.40 | 2.60 | SF | □  |
| YGR286C | BIO2   | 2.60 | 2.93 | 2.24 | 2.59 | NO | □  |
| YJL197W | UBP12  | 2.39 | 3.16 | 2.12 | 2.56 | NO | □  |
| YJR070C | □      | 2.36 | 2.45 | 2.80 | 2.54 | NO | □  |
| YDL188C | PPH22  | 2.59 | 2.73 | 2.25 | 2.53 | SF | SF |
| YDR014W | RAD61  | 2.65 | 2.82 | 2.06 | 2.51 | SF | SF |
| YDR431W | □      | 2.80 | 2.14 | 2.59 | 2.51 | SF | SF |
| YGL249W | ZIP2   | 2.42 | 3.01 | 2.09 | 2.51 | NO | □  |
| YIL049W | DFG10  | 2.30 | 2.91 | 2.30 | 2.50 | NO | □  |
| YKR013W | PRY2   | 2.93 | 2.15 | 2.39 | 2.49 | NO | □  |
| YMR129W | POM152 | 2.39 | 2.63 | 2.46 | 2.49 | NO | □  |
| YJL141C | YAK1   | 2.35 | 2.32 | 2.79 | 2.48 | SF | SF |
| YIL161W | □      | 2.08 | 2.38 | 2.94 | 2.47 | NO | □  |
| YMR127C | SAS2   | 3.08 | 2.17 | 2.04 | 2.43 | NO | □  |
| YKL174C | □      | 3.05 | 2.11 | 2.13 | 2.43 | NO | □  |
| YML028W | TSA1   | 2.41 | 2.48 | 2.39 | 2.43 | NO | □  |
| YIL017C | VID28  | 2.22 | 3.02 | 2.02 | 2.42 | SF | SL |
| YPL155C | KIP2   | 2.44 | 2.50 | 2.32 | 2.42 | SF | SL |
| YGL121C | GPG1   | 2.57 | 2.26 | 2.41 | 2.41 | SF | SF |
| YNL094W | APP1   | 2.14 | 2.74 | 2.34 | 2.41 | NO | □  |

|         |        |      |      |      |      |    |    |
|---------|--------|------|------|------|------|----|----|
| YIL136W | OM45   | 2.73 | 2.18 | 2.26 | 2.39 | NO | □  |
| YOR314W | □      | 2.07 | 2.55 | 2.49 | 2.37 | NO | □  |
| YNR073C | □      | 2.45 | 2.04 | 2.62 | 2.37 | NO | □  |
| YPL055C | LGE1   | 2.68 | 2.31 | 2.07 | 2.35 | SF | SF |
| YDR521W | □      | 2.37 | 2.44 | 2.24 | 2.35 | NO | □  |
| YFR013W | IOC3   | 2.61 | 2.21 | 2.23 | 2.35 | NO | □  |
| YGL138C | □      | 2.03 | 2.23 | 2.77 | 2.34 | NO | □  |
| YIL101C | XBP1   | 2.12 | 2.18 | 2.69 | 2.33 | NO | □  |
| YGR208W | SER2   | 2.23 | 2.51 | 2.20 | 2.31 | NO | □  |
| YER169W | RPH1   | 2.14 | 2.33 | 2.35 | 2.27 | NO | □  |
| YIL167W | □      | 2.02 | 2.44 | 2.36 | 2.27 | NO | □  |
| YIR005W | IST3   | 2.35 | 2.37 | 2.06 | 2.26 | NO | □  |
| YLR433C | CNA1   | 2.02 | 2.35 | 2.40 | 2.26 | NO | □  |
| YIL135C | VHS2   | 2.20 | 2.26 | 2.28 | 2.25 | NO | □  |
| YDR318W | MCM21  | 2.07 | 2.21 | 2.44 | 2.24 | SF | SF |
| YIL124W | AYR1   | 2.01 | 2.26 | 2.44 | 2.24 | NO | □  |
| YMR297W | PRC1   | 2.21 | 2.21 | 2.30 | 2.24 | NO | □  |
| YPL064C | CWC27  | 2.24 | 2.41 | 2.04 | 2.23 | NO | □  |
| YIL140W | AXL2   | 2.00 | 2.29 | 2.23 | 2.17 | NO | □  |
| YJR095W | SFC1   | 2.13 | 2.17 | 2.16 | 2.15 | NO | □  |
| YLR085C | ARP6   | 2.09 | 2.16 | 2.18 | 2.15 | SF | SF |
| YBR273C | □      | 2.09 | 2.07 | 2.26 | 2.14 | NO | □  |
| YPL267W | □      | 2.08 | 2.09 | 2.18 | 2.12 | NO | □  |
| YPL260W | □      | 2.02 | 2.11 | 2.15 | 2.10 | NO | □  |
| YKL053W | □      | 2.00 | 2.02 | 2.20 | 2.07 | NO | □  |
| YPL152W | RRD2   | 2.13 | 1.92 | 3.16 | 2.40 | SF | SL |
| YML124C | TUB3   | 2.06 | 1.64 | 1.90 | 1.87 | SF | SF |
| YJR135C | MCM22  | 1.82 | 1.62 | 2.04 | 1.82 | SF | SF |
| YER016W | BIM1   | 1.52 | 3.05 | 0.77 | 1.78 | □  | □  |
| YCL029C | BIK1   | 1.57 | 1.75 | 1.95 | 1.76 | SF | SL |
| YMR138W | CIN4   | 1.87 | 1.67 | 1.72 | 1.75 | SF | □  |
| YDR254W | CHL4   | 1.38 | 1.54 | 2.06 | 1.66 | SF | SF |
| YGL168W | HUR1   | 1.76 | 1.79 | 1.42 | 1.66 | SF | SF |
| YOR265W | RBL2   | 1.78 | 1.46 | 1.67 | 1.64 | □  | □  |
| YLR381W | CTF3   | 1.45 | 1.22 | 1.56 | 1.41 | □  | □  |
| YER155C | BEM2   | 1.42 | 1.32 | 1.05 | 1.26 | SF | SF |
| YPR046W | MCM16  | 1.31 | 1.07 | 1.31 | 1.23 | □  | SF |
| YBR107C | IML3   | 0.99 | 1.11 | 1.12 | 1.07 | □  | □  |
| YER007W | PAC2   | □    | □    | □    | □    | □  | □  |
| YOR014W | RTS1   | □    | □    | □    | □    | □  | □  |
| YBR122C | MRPL36 | □    | □    | □    | □    | □  | □  |



---

Notes: In each *cin8* dSLAM experiment, useful UPTAG and DOWNTAG ratios were averaged to obtain a single/double (C/E) hybridization ratio for each YKO. A mean ratio was taken among all three experiments and ranked. Confirmation results were obtained by random spore analysis (on synthetic haploid selection media) or tetrad dissection (on YPD) of heterozygous diploid double YKOs or both: SL-synthetic lethal; SF-synthetic fitness defect or synthetic sick; NO-no synthetic interaction; Blank-not individually tested here. Genes colored in blue are unique to dSLAM; genes in red are unique to SGA; genes identified by both screens are in black.

**Table S6. A summary of the *bim1*□ dSLAM screen.**

| ORF       | Gene Name | single/double |          |            | confirmation |        |
|-----------|-----------|---------------|----------|------------|--------------|--------|
|           |           | ratio 1#      | ratio 2# | mean ratio | RSA          | tetrad |
| YJL030W   | MAD2      | 19.57         | 33.10    | 26.34      | SL           | SL     |
| YDR318W   | MCM21     | 28.19         | 17.77    | 22.98      | SL           | SL     |
| YOR058C   | ASE1      | 22.82         | 20.46    | 21.64      | SL           | SL     |
| YJL124C   | LSM1      | 17.42         | 21.65    | 19.53      | SL           | SL     |
| YDR254W   | CHL4      | 13.66         | 20.88    | 17.27      | SL           | SL     |
| YKR061W   | KTR2      | 19.74         | 13.59    | 16.67      | SL           | SL     |
| YMR078C   | CTF18     | 9.47          | 22.02    | 15.74      | SL           | SL     |
| YHR191C   | CTF8      | 18.38         | 12.78    | 15.58      | SL           | SL     |
| YCL029C   | BIK1      | 13.97         | 15.82    | 14.90      | SL           | SL     |
| YPL018W   | CTF19     | 16.00         | 13.14    | 14.57      | SL           | SL     |
| YDR378C   | LSM6      | 10.20         | 16.29    | 13.25      | SL           | SL     |
| YOR073W   | SGO1      | 13.42         | 12.88    | 13.15      | SL           | SL     |
| YOR195W   | SLK19     | 14.35         | 11.49    | 12.92      | SL           | SL     |
| YCL016C   | DCC1      | 11.79         | 13.76    | 12.78      | SL           | SL     |
| YJR135C   | MCM22     | 10.79         | 14.62    | 12.71      | SL           | SL     |
| YDR150W   | NUM1      | 18.07         | 6.56     | 12.32      | SL           | SL     |
| YGL086W   | MAD1      | 9.58          | 14.75    | 12.17      | SL           | SL     |
| YNL273W   | TOF1      | 11.88         | 10.82    | 11.35      | SL           | SL     |
| YOR026W   | BUB3      | 11.95         | 8.81     | 10.38      | SL           | SL     |
| YML094C-A | □         | 10.08         | 10.62    | 10.35      | SL           | SL     |
| YJL013C   | MAD3      | 9.44          | 9.21     | 9.32       | SL           | SL     |
| YBR107C   | IML3      | 9.93          | 8.65     | 9.29       | SL           | SL     |
| YPL174C   | NIP100    | 10.73         | 7.45     | 9.09       | SL           | SL     |
| YLR200W   | YKE2      | 9.09          | 8.77     | 8.93       | SL           | SL     |
| YML094W   | GIM5      | 7.86          | 9.80     | 8.83       | SL           | SL     |
| YMR299C   | DYN3      | 10.31         | 6.88     | 8.60       | SL           | SL     |
| YPL017C   | □         | 5.38          | 11.61    | 8.50       | SL           | SL     |
| YPR135W   | CTF4      | 7.19          | 8.96     | 8.08       | SL           | SL     |
| YGL216W   | KIP3      | 7.49          | 8.62     | 8.05       | SL           | SL     |
| YLR381W   | CTF3      | 8.75          | 7.28     | 8.02       | SL           | SL     |
| YDR149C   | □         | 8.27          | 7.05     | 7.66       | SL           | SL     |
| YGR078C   | PAC10     | 8.59          | 6.39     | 7.49       | SL           | SL     |
| YOL072W   | THP1      | 7.16          | 7.01     | 7.09       | SL           | SL     |
| YLR410W   | VIP1      | 8.81          | 5.16     | 6.98       | SL           | SL     |
| YPL155C   | KIP2      | 6.66          | 6.94     | 6.80       | SL           | SL     |
| YDR159W   | SAC3      | 9.79          | 3.78     | 6.79       | SL           | SF     |
| YDR014W   | RAD61     | 5.43          | 8.12     | 6.78       | SL           | SL     |
| YOR269W   | PAC1      | 9.30          | 3.83     | 6.56       | SL           | SL     |
| YPR046W   | MCM16     | 6.87          | 5.87     | 6.37       | SL           | SL     |

|           |        |      |      |      |    |    |
|-----------|--------|------|------|------|----|----|
| YMR048W   | CSM3   | 5.35 | 7.07 | 6.21 | SL | SL |
| YPL008W   | CHL1   | 5.21 | 7.03 | 6.12 | SL | SL |
| YLL049W   | □      | 6.41 | 5.83 | 6.12 | SL | SL |
| YGL217C   | □      | 5.44 | 6.73 | 6.08 | SL | SL |
| YAL021C   | CCR4   | 6.12 | 6.00 | 6.06 | SF | SF |
| YKR082W   | NUP133 | 3.00 | 8.90 | 5.95 | SL | SL |
| YKR054C   | DYN1   | 6.03 | 5.64 | 5.83 | SL | SL |
| YNL153C   | GIM3   | 5.76 | 5.90 | 5.83 | SL | SL |
| YCR086W   | CSM1   | 4.10 | 6.63 | 5.36 | SL | SL |
| YCR077C   | PAT1   | 6.79 | 3.69 | 5.24 | SL | SL |
| YDR532C   | KRE28  | 4.41 | 5.95 | 5.18 | SL | SL |
| YAL024C   | LTE1   | 4.40 | 5.40 | 4.90 | SL | SL |
| YDR359C   | VID21  | 7.47 | 1.83 | 4.65 | SL | SL |
| YEL003W   | GIM4   | 5.45 | 3.65 | 4.55 | SL | SL |
| YOL076W   | MDM20  | 5.05 | 3.58 | 4.32 | SF | SL |
| YGR180C   | RNR4   | 2.26 | 6.26 | 4.26 | SL | SL |
| YNL107W   | YAF9   | 4.54 | 3.72 | 4.13 | SL | SL |
| YHR129C   | ARP1   | 5.14 | 2.86 | 4.00 | SL | SL |
| YOL012C   | HTZ1   | 4.74 | 3.23 | 3.98 | SL | SL |
| YCR079W   | □      | 4.14 | 3.74 | 3.94 | SF | SF |
| YDL116W   | NUP84  | 2.91 | 4.86 | 3.88 | SF | SL |
| YJR053W   | BFA1   | 3.51 | 4.15 | 3.83 | SL | SL |
| YJL148W   | RPA34  | 3.75 | 3.88 | 3.81 | SL | SF |
| YDR488C   | PAC11  | 3.25 | 4.26 | 3.76 | SL | SL |
| YLR268W   | SEC22  | 3.88 | 3.46 | 3.67 | SF | SF |
| YDR439W   | LRS4   | 3.27 | 4.05 | 3.66 | SL | SL |
| YDR432W   | NPL3   | 3.84 | 3.39 | 3.62 | SL | SL |
| YGR188C   | BUB1   | 2.22 | 4.94 | 3.58 | SL | SL |
| YML041C   | □      | 3.61 | 3.14 | 3.37 | SF | SF |
| YIL153W   | RRD1   | 2.09 | 4.56 | 3.33 | SF | SL |
| YNL171C   | □      | 3.12 | 3.44 | 3.28 | SL | □  |
| YML112W   | CTK3   | 2.14 | 4.39 | 3.26 | SF | SL |
| YOR295W   | UAF30  | 4.52 | 1.77 | 3.14 | SF | SL |
| YPL253C   | VIK1   | 3.39 | 2.82 | 3.11 | SF | SL |
| YDR156W   | RPA14  | 3.07 | 3.08 | 3.07 | SF | SF |
| YER164W   | CHD1   | 3.01 | 3.11 | 3.06 | SF | SF |
| YCR024C   | □      | 4.37 | 1.58 | 2.97 | SF | SF |
| YER014C-A | BUD25  | 3.22 | 2.68 | 2.95 | SF | SL |
| YMR055C   | BUB2   | 3.05 | 2.73 | 2.89 | SL | SL |
| YNR051C   | BRE5   | 2.25 | 3.48 | 2.86 | SF | SL |
| YOL054W   | □      | 3.73 | 1.96 | 2.85 | SF | SF |
| YEL018W   | □      | 3.61 | 2.06 | 2.83 | SF | SF |
| YLR254C   | □      | 2.72 | 2.80 | 2.76 | SL | SL |

|         |        |      |      |      |    |    |
|---------|--------|------|------|------|----|----|
| YGL003C | CDH1   | 2.63 | 2.87 | 2.75 | SL | SL |
| YBR284W | □      | 3.82 | 1.66 | 2.74 | SF | SF |
| YIR005W | IST3   | 2.47 | 2.99 | 2.73 | SF | SF |
| YDR424C | DYN2   | 2.37 | 2.99 | 2.68 | SL | SL |
| YJL179W | PFD1   | 2.77 | 2.55 | 2.66 | SL | SL |
| YFR030W | MET10  | 1.97 | 3.32 | 2.64 | SF | SL |
| YGR270W | YTA7   | 1.59 | 3.70 | 2.64 | SF | SL |
| YDR383C | NKP1   | 2.60 | 2.57 | 2.58 | SF | SF |
| YNL054W | VAC7   | 2.85 | 2.27 | 2.56 | SF | SL |
| YMR293C | □      | 2.58 | 2.52 | 2.55 | SL | SL |
| YJR060W | CBF1   | 3.19 | 1.82 | 2.51 | SL | SL |
| YPR023C | EEF3   | 2.89 | 2.10 | 2.49 | SF | SF |
| YBR231C | AOR1   | 2.32 | 2.63 | 2.47 | SF | SL |
| YER177W | BMH1   | 1.56 | 3.38 | 2.47 | SF | SF |
| YDR334W | SWR1   | 2.64 | 2.30 | 2.47 | SF | SF |
| YKL048C | ELM1   | 2.56 | 2.32 | 2.44 | SF | SF |
| YDL234C | GYP7   | 1.84 | 2.94 | 2.39 | SF | SF |
| YML016C | PPZ1   | 2.05 | 2.66 | 2.36 | SF | SF |
| YKL205W | LOS1   | 1.62 | 3.07 | 2.34 | SF | SF |
| YPL205C | □      | 1.60 | 3.04 | 2.32 | SF | SL |
| YMR294W | JNM1   | 2.12 | 2.49 | 2.30 | SF | SF |
| YOR349W | CIN1   | 2.10 | 2.51 | 2.30 | SF | SF |
| YCR064C | □      | 2.12 | 2.45 | 2.28 | SF | SL |
| YNL225C | CNM67  | 2.33 | 2.14 | 2.23 | SL | SL |
| YKR010C | TOF2   | 2.25 | 2.13 | 2.19 | SF | SF |
| YML103C | NUP188 | 1.77 | 2.60 | 2.18 | SF | SF |
| YPL152W | RRD2   | 2.43 | 1.84 | 2.13 | SF | SL |
| YJL006C | CTK2   | 2.16 | 2.11 | 2.13 | SF | SL |
| YDR360W | □      | 2.54 | 1.68 | 2.11 | SF | SF |
| YNL248C | RPA49  | 1.76 | 2.44 | 2.10 | SF | SF |
| YAR014C | BUD14  | 1.67 | 2.46 | 2.06 | SF | SF |
| YHR013C | ARD1   | 1.90 | 2.17 | 2.03 | SF | SF |
| YIL017C | VID28  | 2.18 | 1.88 | 2.03 | SF | SL |
| YDR174W | HMO1   | 2.48 | 1.59 | 2.03 | SF | SF |
| YNR052C | POP2   | 2.44 | 1.52 | 1.98 | SF | SF |
| YLR079W | SIC1   | 2.20 | 1.74 | 1.97 | SF | SF |
| YLR085C | ARP6   | 2.00 | 1.90 | 1.95 | SF | SL |
| YNL201C | □      | 2.13 | 1.76 | 1.95 | SF | SF |
| YNL136W | □      | 2.22 | 1.65 | 1.93 | SF | SF |
| YGL240W | DOC1   | 2.12 | 1.74 | 1.93 | SF | SL |
| YNL307C | MCK1   | 1.97 | 1.73 | 1.85 | SF | SF |
| YIL040W | □      | 1.78 | 1.90 | 1.84 | SF | SL |
| YIL103W | □      | 1.66 | 1.98 | 1.82 | SF | SF |

|         |        |      |      |      |    |    |
|---------|--------|------|------|------|----|----|
| YKL139W | CTK1   | 2.03 | 1.55 | 1.79 | SF | SL |
| YER095W | RAD51  | 2.03 | 1.53 | 1.78 | SF | SF |
| YDR183W | PLP1   | 1.93 | 1.55 | 1.74 | SF | SF |
| YCL061C | MRC1   | 1.51 | 1.96 | 1.74 | SF | SL |
| YLR399C | BDF1   | 1.86 | 1.56 | 1.71 | SF | SL |
| YCL060C | □      | 1.37 | 6.59 | 3.98 | SF | SL |
| YJR063W | RPA12  | 3.63 | □    | 3.63 | SF | SF |
| YLR315W | NKP2   | 1.24 | 3.14 | 2.19 | SF | SL |
| YIL036W | CST6   | 0.74 | 3.64 | 2.19 | SF | SF |
| YPL200W | CSM4   | 1.09 | 3.13 | 2.11 | SF | SF |
| YGR200C | ELP2   | 3.01 | 1.12 | 2.06 | SF | SF |
| YFR036W | CDC26  | 1.13 | 2.90 | 2.01 | SF | SL |
| YDR289C | RTT103 | 2.57 | 1.41 | 1.99 | SF | SF |
| YOL081W | IRA2   | 0.76 | 3.08 | 1.92 | NO | SF |
| YLR204W | QRI5   | 2.51 | 1.16 | 1.83 | SF | SF |
| YLR182W | SWI6   | 2.64 | 1.00 | 1.82 | SF | SL |
| YNL298W | CLA4   | 2.16 | 1.46 | 1.81 | SF | SF |
| YPL241C | CIN2   | 2.06 | 1.48 | 1.77 | SF | SF |
| YDR315C | IPK1   | 2.28 | 1.20 | 1.74 | SF | SL |
| YFR019W | FAB1   | 2.49 | 0.98 | 1.74 | SF | SL |
| YOL004W | SIN3   | 2.47 | 0.94 | 1.70 | SF | SF |
| YMR267W | PPA2   | 1.92 | 1.47 | 1.70 | SF | SL |
| YEL061C | CIN8   | 0.70 | 2.69 | 1.70 | NO | SF |
| YNL170W | □      | 1.98 | 1.39 | 1.68 | SF | SF |
| YMR073C | □      | 2.12 | 1.14 | 1.63 | SF | SF |
| YNL330C | RPD3   | 1.86 | 1.34 | 1.60 | SF | SL |
| YGL252C | RTG2   | 1.07 | 2.06 | 1.56 | NO | SF |
| YAL010C | MDM10  | 0.93 | 2.15 | 1.54 | NO | SF |
| YJR104C | SOD1   | 1.03 | 2.04 | 1.53 | SF | SL |
| YOR141C | ARP8   | 2.06 | 0.96 | 1.51 | SF | SF |
| YDL040C | NAT1   | 1.63 | 1.29 | 1.46 | SF | SL |
| YPR120C | CLB5   | 1.68 | 1.17 | 1.42 | SF | SF |
| YDL115C | IWR1   | 0.72 | 2.05 | 1.39 | SF | SL |
| YNL296W | KRE25  | 2.14 | 0.24 | 1.19 | SF | SF |
| YML124C | TUB3   | 1.49 | 1.90 | 1.70 | NO | SF |
| YDL117W | CYK3   | 1.25 | 1.48 | 1.36 | □  | □  |
| YLR190W | MMR1   | 1.28 | 1.35 | 1.32 | □  | □  |
| YNL106C | INP52  | 1.11 | 1.27 | 1.19 | NO | NO |
| YPR164W | MMS1   | 1.46 | 1.18 | 1.32 | □  | □  |
| YDL225W | SHS1   | 1.38 | 1.16 | 1.27 | □  | □  |
| YGL163C | RAD54  | 1.15 | 1.13 | 1.14 | NO | SF |
| YMR224C | MRE11  | 0.73 | 1.10 | 0.91 | NO | SF |
| YGR229C | SMI1   | 1.11 | 1.10 | 1.11 | □  | □  |

|         |       |      |      |      |    |    |
|---------|-------|------|------|------|----|----|
| YML032C | RAD52 | 1.12 | 1.08 | 1.10 | NO | SF |
| YGR228W | □     | 0.98 | 1.07 | 1.02 | □  | □  |
| YLR386W | VAC14 | 1.02 | 1.04 | 1.03 | NO | SF |
| YLR373C | VID22 | 1.23 | 1.01 | 1.12 | NO | SF |
| YAL013W | DEP1  | 0.87 | 0.97 | 0.92 | □  | □  |
| YNL215W | IES2  | 0.83 | 0.91 | 0.87 | □  | □  |
| YDR162C | NBP2  | 1.18 | 0.88 | 1.03 | NO | NO |
| YMR263W | SAP30 | 0.94 | 0.82 | 0.88 | NO | SF |
| YBR095C | RXT2  | 1.33 | 0.74 | 1.04 | □  | □  |
| YGL173C | KEM1  | 1.15 | 0.69 | 0.92 | NO | SF |
| YBR200W | BEM1  | 0.74 | 0.69 | 0.71 | NO | SF |
| YGL211W | □     | 0.87 | 0.67 | 0.77 | □  | □  |
| YLR262C | YPT6  | 0.78 | 0.60 | 0.69 | □  | □  |
| YMR312W | ELP6  | 0.82 | 0.59 | 0.70 | NO | SF |
| YKL110C | KTI12 | 1.30 | 0.59 | 0.94 | □  | □  |
| YPL086C | ELP3  | 1.07 | 0.53 | 0.80 | NO | SF |
| YNL097C | PHO23 | 0.69 | 0.49 | 0.59 | NO | SF |
| YLR384C | IKI3  | 1.01 | 0.38 | 0.69 | □  | □  |
| YHR111W | UBA4  | 1.13 | 0.37 | 0.75 | □  | □  |
| YPL102C | KRE24 | 1.04 | 0.35 | 0.70 | □  | □  |
| YDR335W | MSN5  | 1.26 | □    | 1.26 | □  | □  |
| YPL101W | ELP4  | 0.80 | □    | 0.80 | NO | SF |

Notes: In each *bim1*□ dSLAM experiment, useful UPTAG and DOWNTAG ratios were averaged to obtain a single/double (C/E) hybridization ratio for each YKO. A mean ratio was taken between both experiments and ranked. For YKOs with only one data, that data was used without averaging. Confirmation results were obtained by random spore analysis (on synthetic haploid selection media) or tetrad dissection (on YPD) of heterozygous diploid double YKOs or both: SL-synthetic lethal; SF-synthetic fitness defect or synthetic sick; NO-no synthetic interaction; blank-not individually tested here. Genes colored in blue are unique to dSLAM; genes colored in red are unique to SGA; genes identified by both screens are in black.

**Table S7. *sgs1* synthetic interactions identified by dSLAM**

| Gene Name | C/ E mean ratio <sup>a</sup> | RSA confirmation | SGA              | Haploid SLAM | Biological process          |
|-----------|------------------------------|------------------|------------------|--------------|-----------------------------|
| SLX4      | 39.33                        | SL               | Yes <sup>b</sup> | Yes          | DNA replication             |
| YBR100W   | 28.81                        | SL               | Yes              | Yes          | overlaps with MMS4 ORF      |
| MMS4      | 24.26                        | SL               | Yes              | Yes          | DNA repair                  |
| YBR099C   | 17.63                        | SL               | Yes              | Yes          | overlaps with MMS4 ORF      |
| SLX1      | 11.49                        | SL               | Yes              | Yes          | DNA repair                  |
| RTT107    | 9.76                         | SF               | Yes              | Yes          | DNA transposition           |
| WSS1      | 5.49                         | SL               | Yes              | Yes          | protein sumoylation         |
| RAD27     | 2.38                         | SL               | Yes              | Yes          | DNA replication and repair  |
| HPR5      | 8.03                         | SL               | Yes              | No           | DNA repair                  |
| RRM3      | 5.72                         | SF               | Yes              | No           | DNA replication             |
| MUS81     | 4.67                         | SL               | Yes              | No           | DNA repair                  |
| CSM3      | 3.95                         | SF               | Yes              | No           | meiosis                     |
| POL32     | 2.85                         | SF               | Yes              | No           | DNA replication and repair  |
| SOD1      | 2.63                         | SL               | Yes              | No           | superoxide metabolism       |
| ASF1      | 2.53                         | SF               | Yes              | No           | DNA damage response         |
| ESC2      | 2.42                         | SF               | Yes              | No           | chromatin silencing         |
| SAE2      | 2.10                         | SL               | Yes              | No           | DNA damage response         |
| YBR094W   | 2.02                         | SF               | Yes              | No           | unknown                     |
| RAD50     | 1.93                         | SF               | Yes              | No           | DNA repair                  |
| MGS1      | 2.26                         | SF               | Yes              | No           | DNA replication             |
| SLX8      | 14.30                        | SL               | No <sup>c</sup>  | Yes          | DNA recombination           |
| RNH35     | 4.55                         | SF               | No               | Yes          | DNA replication             |
| LYS7      | 3.51                         | SL               | No               | Yes          | copper ion transport        |
| RNH203    | 2.45                         | SF               | No               | Yes          | unknown                     |
| HST3      | 2.14                         | SF               | No               | Yes          | chromatin silencing         |
| HEX3      | 10.45                        | SL               | No               | No           | DNA recombination           |
| NUP84     | 6.33                         | SF               | No               | No           | Nuclear transportation      |
| NUP133    | 5.56                         | SF               | No               | No           | Nuclear transportation      |
| DIA2      | 4.61                         | SF               | No               | No           | invasive growth             |
| RAD26     | 4.44                         | SF               | No               | No           | DNA repair                  |
| YDL162C   | 4.24                         | SF               | No               | No           | unknown                     |
| VID22     | 3.72                         | SF               | No               | No           | vacuolar protein catabolism |
| MRC1      | 3.65                         | SF               | No               | No           | DNA replication checkpoint  |
| TOF1      | 3.36                         | SF               | No               | No           | DNA replication checkpoint  |
| POP2      | 3.03                         | SF               | No               | No           | transcriptional regulation  |
| CTF4      | 2.99                         | SF               | No               | No           | DNA repair                  |
| ASC1      | 2.89                         | SF               | No               | No           | unknown                     |
| GAS5      | 2.88                         | SF               | No               | No           | unknown                     |
| YCL060C   | 2.80                         | SF               | No               | No           | unknown                     |

|         |      |    |    |    |                               |
|---------|------|----|----|----|-------------------------------|
| FAB1    | 2.72 | SF | No | No | phospholipid metabolism       |
| SKT5    | 2.69 | SF | No | No | cell wall chitin biosynthesis |
| YDR520C | 2.63 | SF | No | No | unknown                       |
| NPL3    | 2.62 | SF | No | No | mRNA-nucleus export           |
| PAT1    | 2.43 | SF | No | No | unknown                       |
| SPH1    | 2.40 | SF | No | No | actin filament organization   |
| UBC4    | 2.39 | SF | No | No | protein ubiquitination        |
| YOR024W | 2.38 | SF | No | No | unknown                       |
| DUN1    | 2.38 | SF | No | No | DNA damage response           |
| EXO1    | 2.31 | SF | No | No | DNA mismatch repair           |
| CSM1    | 2.29 | SF | No | No | chromosome segregation        |
| NSG2    | 2.19 | SF | No | No | unknown                       |
| NUP170  | 2.18 | SF | No | No | Nuclear transportation        |
| UME6    | 2.09 | SF | No | No | entry into meiosis            |
| VIP1    | 2.08 | SF | No | No | actin cytoskeleton            |
| MTF1    | 2.03 | SF | No | No | transcriptional regulation    |
| PIF1    | 2.01 | SF | No | No | DNA recombination             |
| MMS1    | 1.99 | SF | No | No | DNA repair                    |
| LSM1    | 1.98 | SF | No | No | mRNA catabolism               |
| UAF30   | 1.88 | SF | No | No | transcriptional regulation    |
| XRS2    | 1.81 | SF | No | No | DNA repair                    |

Notes: <sup>a</sup> In each *sgs1* dSLAM experiment, useful UPTAG and DOWNTAG ratios were averaged to obtain a single/double (C/E) hybridization ratio for each YKO. Each of the interactions listed was individually re-tested by random spore analysis (on synthetic haploid selection media): SL-synthetic lethal; SF-synthetic fitness defect or synthetic sick. <sup>b</sup> “Yes” means that synthetic interaction was identified by that method; <sup>c</sup> “No” means that synthetic interaction was not identified by that method.



**Table S8. Potential *cdc101-1* synthetic interactions**

| Gene Name<br>□ | control/experiment mean ratio <sup>a</sup> |                         |                           |                         |
|----------------|--|-------------------------|---------------------------|-------------------------|
|                | <i>cdc102-1</i><br>(30°C)                  | <i>CDC102</i><br>(30°C) | <i>cdc102-1</i><br>(28°C) | <i>CDC102</i><br>(28°C) |
| MDM39          | 12.17                                      | 1.27                    | 7.09                      | 1.51                    |
| NHX1           | 11.40                                      | 1.20                    | 6.27                      | 1.25                    |
| RMD7           | 9.06                                       | 1.41                    | 12.00                     | 1.40                    |
| NBP2           | 8.09                                       | 1.26                    | 2.88                      | 1.08                    |
| YDR455C        | 7.50                                       | 0.96                    | 2.50                      | 1.19                    |
| PTC1           | 6.64                                       | 1.32                    | 2.58                      | 0.81                    |
| ARR4           | 6.16                                       | 1.22                    | 2.68                      | 1.41                    |
| APL5           | 5.96                                       | 1.33                    | 1.98                      | 1.59                    |
| ARL1           | 5.94                                       | 1.32                    | 2.44                      | 0.84                    |
| VPS25          | 5.58                                       | 1.50                    | 2.75                      | 1.17                    |
| YGL072C        | 5.27                                       | 1.26                    | 2.12                      | 1.30                    |
| NKP2           | 5.14                                       | 0.99                    | 1.73                      | 0.74                    |
| YCR051W        | 5.01                                       | 0.48                    | 1.70                      | 1.27                    |
| VPS36          | 4.96                                       | 1.03                    | 2.18                      | 1.08                    |
| RIM20          | 4.96                                       | 1.46                    | 2.27                      | 1.40                    |
| YER084W        | 4.87                                       | 1.40                    | 2.46                      | 1.18                    |
| SPF1           | 4.86                                       | 1.28                    | 2.62                      | 1.03                    |
| MON2           | 4.82                                       | 1.32                    | 2.41                      | 0.78                    |
| APM1           | 4.75                                       | 1.35                    | 2.62                      | 1.71                    |
| HTZ1           | 4.66                                       | 1.50                    | 3.02                      | 1.49                    |
| APL6           | 4.34                                       | 1.31                    | 2.04                      | 1.58                    |
| VTC4           | 4.18                                       | 1.24                    | 1.67                      | 1.18                    |
| CUE1           | 4.15                                       | 1.67                    | 1.95                      | 1.38                    |
| COG8           | 4.08                                       | 1.58                    | 2.68                      | 1.19                    |
| ERG4           | 4.07                                       | 1.05                    | 1.93                      | 1.09                    |
| YFR024C        | 4.05                                       | 1.62                    | 2.31                      | 1.27                    |
| APS3           | 3.98                                       | 1.15                    | 2.20                      | 1.47                    |
| WHI2           | 3.89                                       | 1.56                    | 1.67                      | 0.80                    |
| PMP3           | 3.83                                       | 0.81                    | 3.84                      | 0.66                    |
| SMI1           | 3.82                                       | 0.61                    | 2.07                      | 1.24                    |
| GYP1           | 3.81                                       | 1.54                    | 2.27                      | 1.43                    |
| SAC7           | 3.80                                       | 1.03                    | 2.03                      | 1.18                    |
| API2           | 3.73                                       | 1.32                    | 2.04                      | 1.20                    |
| IMP2'          | 3.71                                       | 0.99                    | 2.07                      | 0.99                    |
| YER119C-A      | 3.68                                       | 1.37                    | 2.05                      | 1.27                    |
| YKL136W        | 3.66                                       | 1.37                    | 2.26                      | 1.30                    |
| LSM1           | 3.63                                       | 1.44                    | 2.60                      | 1.20                    |
| APM3           | 3.59                                       | 1.29                    | 1.72                      | 1.21                    |
| UBP14          | 3.46                                       | 1.24                    | 2.30                      | 1.09                    |

|         |      |      |      |      |
|---------|------|------|------|------|
| COG6    | 3.45 | 1.34 | 1.97 | 1.17 |
| PFD1    | 3.41 | 1.26 | 2.81 | 1.24 |
| RIM21   | 3.35 | 1.39 | 2.26 | 1.49 |
| VPS71   | 3.28 | 1.13 | 2.46 | 1.01 |
| SYS1    | 3.27 | 1.66 | 2.21 | 1.32 |
| YLR057W | 3.24 | 1.22 | 1.91 | 1.45 |
| VAM10   | 3.23 | 1.26 | 1.85 | 1.19 |
| SCS2    | 3.17 | 1.50 | 2.28 | 1.32 |
| YDR334W | 3.15 | 1.23 | 2.11 | 1.65 |
| OPI3    | 3.13 | 1.59 | 2.43 | 0.83 |
| CTK3    | 3.12 | 0.86 | 2.41 | 0.68 |
| GEF1    | 3.11 | 1.11 | 1.86 | 1.38 |
| APL2    | 3.11 | 1.19 | 1.88 | 1.05 |
| CKB2    | 3.09 | 1.30 | 2.24 | 0.91 |
| LTE1    | 3.05 | 1.02 | 2.69 | 1.20 |
| RER1    | 3.02 | 0.81 | 2.00 | 1.16 |
| RIM101  | 3.02 | 1.21 | 1.66 | 1.22 |
| GSG1    | 3.00 | 1.80 | 2.01 | 1.26 |
| OST3    | 2.99 | 1.13 | 1.83 | 1.54 |
| INP51   | 2.97 | 1.09 | 2.06 | 1.15 |
| RPA14   | 2.96 | 0.19 | 1.61 | 1.33 |
| YLR374C | 2.93 | 1.50 | 2.23 | 1.14 |
| SEC72   | 2.90 | 0.88 | 1.78 | 0.89 |
| KTR3    | 2.89 | 1.33 | 1.54 | 1.14 |
| YBT1    | 2.88 | 1.56 | 3.21 | 1.32 |
| MCK1    | 2.86 | 1.31 | 1.95 | 1.22 |
| APS1    | 2.86 | 1.27 | 1.81 | 1.48 |
| SAC1    | 2.85 | 1.24 | 1.87 | 0.86 |
| YKL207W | 2.85 | 1.24 | 2.54 | 1.31 |
| SNF8    | 2.81 | 0.89 | 1.83 | 0.87 |
| KHA1    | 2.81 | 1.17 | 1.71 | 1.04 |
| ALG8    | 2.79 | 1.08 | 1.50 | 1.26 |
| VPS67   | 2.77 | 0.94 | 2.03 | 1.06 |
| YGR122W | 2.67 | 1.65 | 2.46 | 1.76 |
| SWC1    | 2.66 | 1.09 | 3.20 | 1.85 |
| RIM13   | 2.65 | 1.39 | 1.85 | 1.10 |
| RIC1    | 2.65 | 1.09 | 1.84 | 0.80 |
| LGE1    | 2.64 | 1.02 | 2.12 | 1.42 |
| YGL046W | 2.62 | 0.82 | 1.51 | 1.30 |
| YLR294C | 2.60 | 0.74 | 1.84 | 0.91 |
| YDL072C | 2.58 | 1.18 | 1.60 | 1.31 |

Notes: <sup>a</sup> Each of the heterozygous double YKO pools containing the indicated plasmid-borne alleles of *CDC102* were sporulated and selected to generated a *MATa* haploid pool of both single and double deletion mutants (control) and a *MATa* haploid double deletion pool (experiment). The abundance of each *kanMX* YKO mutant in both pools was compared by microarray hybridization of the TAGs. A high C/E ratio indicates synthetic interactions between the corresponding YKO and the *CDC102* allele.

**Table S9. PCR primers for amplifying the gene-disruption cassettes.**

| Cassette                | Forward primer sequence              | Reverse primer sequence              |
|-------------------------|--------------------------------------|--------------------------------------|
| <i>bim1</i> □::URA3     | 5' ATTAGACCAACTATAC<br>GCCGAAGG3'    | 5' AACGTGCCCTTCTTGTC<br>AGTATCA3'    |
| <i>bub2</i> □::URA3     | 5' CGGAGGGTCTTTATCG<br>GCGACT3'      | 5' TACCGTCTTTAATGTGA<br>TCCTAACC3'   |
| <i>cdc102</i> □::natMX  | 5' CCATTAATCCTTTTTAT<br>GTAGTCATGC3' | 5' AAGGATGATCAGTATG<br>CGAAATATCTG3' |
| <i>cin8</i> □::URA3     | 5' GACTTTGTTGAGAGAG<br>ACTCTATTG3'   | 5' CCGAATGCCTCATACCA<br>ATTGATG3'    |
| <i>dot1</i> □::URA3     | 5' TTGTTTTCCCTTGCGTGC<br>CCAATG3'    | 5' TGTGCATACCGCTACTT<br>CTTCATC3'    |
| <i>dun1</i> □::URA3     | 5' GTTGGTTCAATGCAAA<br>TTATTACATC3'  | 5' AGAGAGATTGAATCTTA<br>AAATTGACG3'  |
| <i>ira2</i> □::URA3     | 5' GACAATTCATTGTATT<br>TCCAGTTGTTT3' | 5' CCGTCTGAGTTCATATA<br>TAATTCAC3'   |
| <i>mad2</i> □::URA3     | 5' TGCATAGTCTAACCTA<br>TAGTATTTAG3'  | 5' GCCGAAGTACCCATCTA<br>ATGAAAC3'    |
| <i>mre11</i> □::URA3    | 5' GGATATCACTCAAGCA<br>GAAAAAATCA3'  | 5' ATGGGCTTTCTGGCCTG<br>ATCAAC3'     |
| <i>rad51</i> □::URA3    | 5' TCGTACTAATTGCTGA<br>CAAACCCG3'    | 5' CTAGAAGGTACTGTGCA<br>TCAGGG3'     |
| <i>rad61</i> □::URA3    | 5' TCTTTGAACTTGTTGC<br>AATTGG3'      | 5' CCTCATGTTGTCTAACTT<br>CCAA3'      |
| <i>rad9</i> □::URA3     | 5' AGCTCTTGAACAACAT<br>ACTCTCAG3'    | 5' GAGATTCATCAAACAGA<br>TTGATCGC3'   |
| <i>sgs1</i> □::URA3     | 5' GGTAATAGCGCGATGA<br>AACAACGTC3'   | 5' GTGCTCAAGGTGAATAT<br>ACTGGTC3'    |
| <i>sml1</i> □::URA3     | 5' CATGATGGTGTAAGTC<br>ACAAAGAC3'    | 5' ACGATACGTGAAGGCAT<br>ATGAA3'      |
| <i>tub1</i> □::URA3     | 5' GTGCTTGCCAACACCC<br>ATATT3'       | 5' GTCCATATTAAATCCGC<br>TGATAG3'     |
| <i>ura3</i> □::URA3     | 5' TCCCGGTTAGTACAGA<br>TGC3'         | 5' GATGTGGTGCTGATTCA<br>GG3'         |
| <i>ykl088w</i> □::natMX | 5' TGATTTTGCCACACCA<br>ATCACACAG3'   | 5' CTGCTAGTGCCTACTCT<br>AGTTGTTG3'   |
| <i>ykr079c</i> □::natMX | 5' GGTATCTGACATTATG<br>GTACAGGTT3'   | 5' GGGCACTGGACAAAAA<br>GATAACATT3'   |
| <i>yor145c</i> □::natMX | 5' TCAAACCTTCGAGCTTG<br>CCATCCTAT3'  | 5' CACAAAGCTTGCTGATC<br>ATTGCCCA3'   |

**Reference**

- S1. Bennett, C. B., Lewis, L. K., Karthikeyan, G., Lobachev, K. S., Jin, Y. H., Sterling, J. F., Snipe, J. R., and Resnick, M. A. (2001). Genes required for ionizing radiation resistance in yeast. *Nat Genet* 29, 426-434.
- S2. Chang, M., Bellaoui, M., Boone, C., and Brown, G. W. (2002). A genome-wide screen for methyl methanesulfonate-sensitive mutants reveals genes required for S phase progression in the presence of DNA damage. *Proc Natl Acad Sci U S A* 99, 16934-16939.